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Contents | Inhalt

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1. Contents

1. Contents	2	6. Function	23
1.1 Contents	3	6.1 General Functional Description	23
1.2 FCC	5	6.2 Test Programs	23
2. Overview and Proper Use	6	6.3 Managing the Test Results	24
2.1 Overview	6	6.4 Function Test	24
2.2 Proper Use	6	6.5 Cleaning	24
2.3 Technical Specifications	7	6.6 User Interrupt	24
2.3.1 Dimensions and Weight	7	6.7 Self-Test for Broken Conductors and Short-Circuiting	24
2.3.2 Connection Data	7	6.8 Calibration	24
2.3.3 Pneumatics	7	6.9 Internal Pressure Sensor	24
2.3.4 Test Media	7	6.10 Pressure Gauge Function	24
2.3.5 Product-specific Data	8	6.11 Cleaning Drying	24
2.3.6 Test Methods	9	6.12 Safety Functions	24
2.3.7 General Data	9	6.13 Software Update	25
3. Safety	10	6.14 Serial Number	25
3.1 Notes Warnings Descriptions	10	7. Operation	26
3.1.1 Manufacturer's ID Label	10	7.1 Switching on the Unit	26
3.2 Installed Safety Equipment	11	7.1.1 Operating Elements	26
3.3 Safety Precautions	11	7.2 Switching the Unit ON and OFF	28
3.4 Operator's Obligations	12	7.2.1 Running a Test	29
3.5 Safety Tests and Inspections Carried Out by the Manufacturer at the Facilities	12	7.2.2 Switching Off the Unit	29
4. General Warnings	13	7.2.3 Switching Off the Unit in the Event of a Power Failure	29
4.1 Hazards	13	7.3 Tests	29
4.2 Operating Personnel and Service and Maintenance Personnel	13	7.3.1 Test Parameters	29
4.3 Installation of Replacement Parts	14	7.3.2 Test Methods	30
4.4 Shutdown Procedure	14	7.3.3 Programming Tests (F2 – Main Menu)	30
5. Installation	15	7.4 Diffusion Test	31
5.1 Equipment Supplied	15	7.5 Bubble Point Test	35
5.2 Receipt of Delivery	15	7.6 Customer-Specific BP	36
5.3 Installation Instructions	16	7.7 Diffusion and Bubble Point Test (Combined Test)	36
5.3.1 Ambient Conditions and Conditions at the Place of Installation	16	7.8 Water Intrusion Test Water Flow Test	37
5.3.2 Transport and Unpacking the Equipment	16	7.9 Pressure Drop Test	38
5.4 Connections	17	7.10 Multipoint Diffusion Test	38
5.4.1 Front Panel and Back Panel	17	7.11 Customer-Specific Multipoint Diffusion Test	39
5.4.2 Left and Right Panels	17	7.12 Volume Measurement	39
5.4.3 Left and Right Panels Sartochek® 4 MultiUnit	18	7.13 Plate Holder Test	40
5.5 Connections with Internal Pressure Sensor	18	7.14 Installed Bag Test	41
5.6 Connections with External Pressure Sensor	20	7.15 Program Database (F4 – Main Menu)	42
5.6.1 Inserting Paper Roll and Ribbon Cartridge	21	7.15.1 Loading Test Program from Internal Memory (F1)	43
5.6.2 Sartochek® 4 MultiUnit	22	7.15.2 Uploading Test Program from SD card (F2)	43
		7.15.3 Saving Backup to SD card (F4)	45
		7.15.4 Loading Backup from SD card (F6)	46
		7.15.5 Deleting PLC Programs (F8)	47
		7.16 Visualize (F3 – Main Menu)	47

7.17	Results Database (F5 – Main Menu)	47	13. Electronic Data Transfer via Network Connection . .	84
7.18	Special Functions (F6 – Main Menu)	49	13.1	Network Concept 84
7.18.1	Function Test MU (F1)	49	13.2	FTP Server Requirements 85
7.18.2	Function Test TU (F2)	50	13.2.1	How to Configure Your FTP Server 85
7.18.3	Cleaning (F3)	50	14. PLC Control	86
7.18.4	Drying (F4)	52	14.1	Digital Inputs 86
7.18.5	Pressure Gauge Function (F5)	52	14.1.1	START Input 86
7.18.6	Modify Own Password (F6)	53	14.1.2	PROG Input 86
7.19	Service (F7 – Main Menu)	53	14.1.3	STOP Input 86
7.19.1	Load Program from TU (F1-Service)	54	14.2	Digital Output 86
7.19.1.1	External reference volume	55	14.2.1	POWER on 86
7.19.1.2	Network Settings	56	14.2.2	TEST RUNNING 86
7.19.1.3	Test Unit Address Allocation	57	14.2.3	ERROR 86
7.19.2	Hardware Configuration (F2-Service)	58	14.2.4	TEST OK 86
7.19.3	Calibration Adjustment (F3-Service)	62	14.2.5	TEST not OK 86
7.19.3.1	Calibration and Adjustment of the Internal Pressure Sensor	63	14.2.6	External Valve 1 86
7.19.3.2	Calibration and Adjustment of the External Pressure Sensor Input	65	14.2.7	External Valve 2 86
7.19.4	User Data Management (F4-Service)	66	14.3	Pin Layout 87
7.20	Using the Barcode Scanners (optional)	67	14.4	Examples 88
7.20	Using the Barcode Scanners (optional)	67	14.5	Outlets 91
8.	Maintenance Servicing	68	15. Index	92
8.1	Maintenance	68		
8.1.1	Exchanging the Fine-wire Fuses	68		
8.1.2	Sartorius Stedim Biotech Service Menu	68		
8.2	Cleaning	68		
9.	Troubleshooting	69		
9.1	Disposal and repair instructions	69		
9.2	Troubleshooting	69		
10.	Emergencies	80		
11.	Information and Instructions on Disposal and Repairs	80		
12.	Appendix	81		
12.1	Detailed Test Description	81		
12.1.1	Diffusion Test (Only valid for the Sartocheck® 4 plus Filter tester)	81		
12.1.2	Bubble Point Test (Only valid for the Sartocheck® 4 plus Filter tester)	82		
12.1.3	Pressure Drop Test	83		
12.1.4	Plate Holder Test	83		

1.2 FCC

Note:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference using one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio | TV technician for help.

Shielded Cables

Connections between the devices and peripherals must be made using shielded cables in order to maintain compliance with FCC radio frequency emission limits.

Modifications

Any modifications made to this device that are not approved by Sartorius Stedim Biotech GmbH may void the authority granted to the user by the FCC to operate this equipment.

This instruction manual is based on screen shots from software version 3.1.3. It is also valid for later software versions.

Translation

When the unit is supplied to countries in the EEA, the operating instructions must be translated into the language or languages of the country of use. If any discrepancies arise in the translated text, the original German operating instructions shall form the basis for clearer interpretation or the manufacturer must be contacted.

Copyright

Reproduction, dissemination, transfer by any means and other use is not permitted without the prior written consent of Sartorius Stedim Biotech GmbH. Offenders will be liable for damages. All rights reserved.

2. Overview and Proper Use

2.1 Overview: Sartochek® 4 plus
Filter tester

2. Overview and Proper Use

Fig. 2-1 | Overview

Picture of Sartochek® 4 plus Filter tester

- 1 Housing
- 2 Printer
- 3 Touch screen
- 4 Keypad
- 5 Indicator lights

Overview: Sartochek® 4 plus Filter tester MultiUnit

Fig. 2-2 | Overview

2.2 Proper Use

Sartochek® 4 plus Filter tester is a filter integrity tester. It was exclusively developed, constructed and built for the industrial and commercial purpose of conducting filter tests in pharmaceutical and biotechnological production and laboratory operations.

The unit is not intended for use in or with explosive atmospheres!

The system is exclusively designed for the purpose defined above. Any other application above or beyond this or any modifications to the unit without the written permission of the manufacturer is not considered proper use.

Sartorius Stedim Biotech GmbH cannot assume liability if Sartocheck® 4 plus test units are subjected to improper use or improper operation. Sartorius Stedim Biotech is not liable for any damages resulting therefrom. The operator bears the sole risk.

Proper use also includes compliance with operating, inspection and maintenance conditions prescribed by the manufacturer. Non-compliance with the warnings in the operating instructions and installation instructions, in particular, constitutes improper use of the equipment. Only EMC-tested control leads and accessories may be used. Units may only be returned in their original boxes or appropriate secure packaging. No liability can be accepted for transport damage.

2.3 Technical Specifications

2.3.1 Dimensions and Weight

Sartocheck® 4 plus Filter tester

Length:	460 mm
Width:	390 mm
Height:	245 mm
Weight:	15 kg
Packing Weight:	2.7 kg

Sartocheck® 4 plus Filter tester MultiUnit

Length:	460 mm
Width:	390 mm
Height:	150 mm
Weight:	13 kg
Packing Weight:	2.7 kg

2.3.2 Connection Data

Sartocheck® 4 plus Filter tester | Sartocheck® 4 plus Filter tester MultiUnit

Power requirements

Operating voltage:	100-240 V/±10%
Number of phases:	1 Ph/N
Neutral conductor:	Loadable N (in the Appendix)
Grounding conductor:	Separate SL (yellow-green) (in the Appendix)
Frequency:	47-63 Hz
Protection rating:	IP 44
Installation:	According to VDE regulations
Permissible voltage fluctuation:	± 10% from normal value
Maximum operating pressure:	10 bar g
Operating pressure:	min. 4 bar g
Maximum power input:	74 Watt
Internal limit pressure:	10 bar g

2.3.3 Pneumatics

Sartocheck® 4 plus Filter tester | Sartocheck® 4 plus Filter tester MultiUnit

Required compressed air	4-9 bar g
Water, oil and particle free	

2.3.4 Test Media

Sartocheck® 4 plus Filter tester | Sartocheck® 4 plus Filter tester MultiUnit

Test gas:	Compressed air, dehydrated and deoiled
Wetting fluid:	Water or as described with the process validation documents for the filter cartridge to be tested

2.3.5 Product-specific Data

Manual input range upstream system volume 0¹⁾ – 150000 ml

Measuring Ranges

- Test pressure: 100–8000 mbar g
- Pressure drop: 1–2000 mbar
- Max. measurable volume: 150 L
- Diffusion: 0.1–3000 ml/min
- Water intrusion: 0.1–999 ml/10 min
- Water flow: 0.01–99.9 ml/min
- Bubble point: 100–8000 mbar

Measuring Accuracy

- Pressure $\pm 0.1\% \mid \pm 9.5$ mbar full scale
- Pressure drop ± 1 mbar
- Net volume measurement $\pm 4\%$
- Diffusion $\pm 5\%$
- Water intrusion $\pm 5\%$
- Bubble point ± 50 mbar

Operating Conditions

- Room temperature $+15^{\circ}\text{C}$ to $+35^{\circ}\text{C}$
- Rel. humidity 10–80%

Touch screen (Sartocheck® 4 plus Filter tester only)

- Size 10.4"
- Color 256 colors

Interface Port for Data Output

- Serial port TU RS232 | RS485
- Serial port MU RS232 (Sartocheck® 4 plus Filter tester only)
- PLC port binary signals 12 pins
- Network Connection RJ45

Selectable languages:

English, German, French, Spanish, Italian

¹⁾ When 0 is input, the upstream system volume is measured. Volume input and volume measurement are not applicable for bag testing.

2.3.6 Test Methods

Test method	Without ext. ref. tank	With ext. ref. tank
Press. drop test	All filter systems up to net volume of 150 L	Not applicable
Diffusion test	Cartridge systems up to net volume of 14 L	Cartridge systems up to net volume of 150 L
Bubble point test	Small systems and cartridge systems up to net volume of 14 L	Small systems and cartridge systems up to net volume of 150 L
Multipoint diffusion test	Cartridge systems up to net volume of 14 L	–
Diffusion and bubble point test (complete test)	Cartridge systems up to net volume of 14 L	Cartridge systems up to net volume of 150 L
Intrusion test, water flow test	Cartridge systems up to net volume of 14 L	–

2.3.7 General Data

Ambient conditions	
Ambient temperature:	+5 up to +35°C
Noise level:	< 70 dB (A)
Humidity:	95% at 25°C

The unit should be set up under ambient conditions that produce no negative effects during operation and processing, e.g. from contaminated airflows, metallic dust or fluids from leaks, condensation, aerosols or widely changing temperature conditions.

3. Safety

3.1 Notes | Warnings |
Descriptions

Operating instructions are binding;
designated by a "book".


Warnings
are designated by a "Hazard" triangle.

Electrical Hazards
are designated by the symbol to the left.

Notes
are designated by a "Hand".

Grounding Conductor Connection
is designated by this symbol at the connection sites.

Mandatory Instructions
for handling electrical plug connections.

3.1.1 Manufacturer's ID Label

Fig. 3-1 | Manufacturer's ID Label Sartocheck® 4 plus Filter tester

Fig. 3-2 | Manufacturer's ID Label Sartocheck® 4 Filter tester MultiUnit

The information given in these operating instructions applies to Sartocheck® 4 plus Filter tester, Order no. 26288---FT and Sartocheck® 4 Filter tester MultiUnit, Order no. 26288---FM.

The manufacturer's label is affixed to the rear panel of the unit.

Important: if you have any questions

Always have the following information ready:

- Model
- Version
- Serial number

Having the right information readily available will ensure correct and quick processing.

3.2 Installed Safety Equipment

The installed safety equipment must be inspected at regular test intervals (**d** = daily, **w** = weekly, **m** = monthly). The test methods used are: **S** = sight inspection, **F** = functional testing.

The unit has the following safety equipment installed:

Test	
Interval	Method
m	S/F

Line Disconnecter

The unit is disconnected from the mains voltage by switching off the unit or when the unit is unplugged from the wall outlet (mains supply).

Test	
Interval	Method
m	S

Housing

All unit-specific components are shielded against external influences by a plastic housing. It is strictly forbidden to disengage or dismantle the safety equipment. The function of this safety equipment must be inspected and tested according to the instructions given in the Chapter "Maintenance".

The operator is responsible for instructing the operating personnel and service and maintenance personnel about the meaning and the function of the safety equipment.

These operating instructions are an integral part of the unit and must be available to operating personnel and service and maintenance personnel at all times.

3.3 Safety Precautions (to be carried out by the operator)

The safety precautions and warnings contained herein must be observed.

The operator must

- instruct their operating and servicing personnel in the use of the filter cartridge housing's safety devices and
- ensure that the safety precautions are complied with.
- prevent unauthorized persons (not operating and servicing personnel) from accessing the filter cartridge housing danger zone.

These operating instructions must be kept in a safe place for future reference. The frequency of inspections and checks must be observed.

The work described in these operating instructions is explained in such a way that it can be understood by trained workers and qualified technicians.

Definitions Adapted from EN 60204-1:

Trained Worker

A person familiarized with and, if applicable, trained in their particular duties by a qualified technician and instructed about the potential dangers in the event of improper conduct and on the necessary safety devices and safety precautions.

Qualified Technician

A person who, by virtue of their technical training, know-how, experience and knowledge of relevant standards, is able to assess the work assigned to them and identify potential dangers.

3.4 Operator's Obligations

In the European Economic Area (EEA), the laws to nationally implement Council Directive (89/391/EEC), as well as the corresponding individual Directives and specifically including Council Directive (89/655/EEC) concerning the minimum safety and health requirements for the use of work equipment by workers during labor, must be observed and complied with in their respectively valid versions.

In Germany, the Operational Safety Regulation dated October 2002 must be observed (regulation transposing the above-mentioned Directives into national law). The operator must obtain a local operating permit and follow the regulations and restrictions set forth therein.

Additionally, the operator must comply with the local rules and regulations for:

- The safety of the personnel (accident prevention regulations)
- The safety of the work equipment (protective gear and equipment and maintenance)
- Product disposal (Waste Management Act)
- Material disposal (Waste Management Act)
- Cleaning (cleaning agents and disposal) and the environmental rules and regulations.

3.5 Safety Tests and Inspections Carried Out by the Manufacturer at the Facilities

1. Measurement of airborne noise
2. Testing and inspection according to DIN EN 61010-1

4. General Warnings

4.1 Hazards

4. General Warnings

The safety systems and safety notes described in these operating instructions should be complied with accordingly.

Keep the work area of the Bag tester free of objects during operation so that unimpeded access is possible at all times.

Before any installation, servicing or repair work is carried out, make sure you are aware of the electrical hazards!

Unplug the equipment from the wall outlet (mains supply) before opening it!

Observe the required operating conditions, especially with regard to the temperature of the medium. It is transmitted to the components. Touching components or contact with the medium can result in burns or scalding.

In order for the Sartocheck® 4 plus Filter tester to be EMC compliant, the original Ethernet cable must be used.

4.2 Operating Personnel and Service and Maintenance Personnel

Operating personnel and service and maintenance personnel are persons who are responsible for the transport, installation, operation, equipping, maintenance and cleaning as well as troubleshooting on the unit.

1. The unit may only be operated by trained and authorized persons in accordance with operating instructions that have been issued by the operator.
2. The responsibilities for operation of the unit must be clearly set down and observed so that no ambiguities arise with respect to safety.
3. The safety precautions and warnings given in the operating instructions must be complied with when carrying out any work (operation, maintenance, repairs etc.).
4. The operator shall never work in any way that could impair the safety of the unit or him | herself or others.
5. It is the responsibility of the operator to ensure that only authorized persons work on the unit.
6. It is the responsibility of the operator to report immediately any changes occurring on or to the unit to the operating company.
7. It is the responsibility of the operator to ensure that the unit is only operated in perfect working condition.

4.3 Installation of

Please note that replacement parts and Replacement Parts accessories that not have been supplied by Sartorius Stedim Biotech GmbH have not been tested and approved by Sartorius Stedim Biotech GmbH. The installation of and | or the use of such products can therefore have a negative effect on the design properties of the unit as constructed.

Sartorius Stedim Biotech GmbH is not liable for any damages resulting from the use of non-original parts or non-original accessory parts.

Standard parts can be purchased from specialized dealers.

4.4 Shutdown Procedure

The following shutdown procedure must be followed before any cleaning, servicing or repair work is carried out.

1. Line disconnect
Turn off the main switch by pulling out the AC adapter of the mains supply.
2. Compressed air supply
Turn off the compressed air supply and vent system.

Personnel may suffer physical injuries or fatalities if this procedure is not followed!

5. Installation

5.1 Equipment Supplied

5. Installation

Fig. 5-1 | Sartocheck® 4 plus Filter tester

The equipment supplied with the Sartocheck® 4 plus Filter tester includes:

1. Sartocheck® 4 plus Filter tester, Order no. 26288---FT
 2. Operating instructions and technical documents
 3. Accessories (not shown in the picture)
- Tubing "Inlet" with compressed air filter (I.D. = 6 mm) Order no.: 18104
 - Tubing "Outlet" (I.D. = 4 mm) Order no.: 18103
 - Spare ink ribbon cartridge Order no.: 6982141
 - Spare roll of printer paper Order no.: 6982142
 - Country-specific power cord

The equipment supplied with the Sartocheck® 4 Filter tester MultiUnit includes:

1. Sartocheck® 4 Filter tester MultiUnit, Order no. 26288---FM
 2. Operating instructions and technical documents
 3. Accessories (not shown in the picture)
- Tubing "Inlet" with compressed air filter (I.D. = 6 mm) Order no.: 18104
 - Tubing "Outlet" (I.D. = 4 mm) Order no.: 18103
 - Country-specific power cord
 - Connecting cable Sartocheck® 4 plus Filter tester MultiUnit and Sartocheck® 4 plus Bag tester MultiUnit, Order no.: 1ZE---0010

5.2 Receipt of Delivery

Incoming inspection:

- Check the shipment for damage!
- Check the shipment against the delivery note to ensure that it is complete!

Complaints

If the shipment was damaged during transit:

- Immediately contact the carrier!
- Keep the packaging material (in case the carrier wants to inspect it or the shipment has to be returned to us)

If the shipment is not complete:

- Immediately contact Sartorius Stedim Biotech GmbH or the Customer Service Department in your country!

If the delivery has other defects:

- Immediately contact Sartorius Stedim Biotech GmbH or the Customer Service Department in your country!

Packaging for Return Shipment

Use the original packaging and original packaging material as far as possible.

If the original packaging and original packaging material are no longer available:

- Commission a packaging company with specialized technicians.
- Pack the unit in a shipping container in such a way that mechanical damage is prevented.

If you have any questions about packaging and transport securing devices, please contact Sartorius Stedim Biotech GmbH.

Check the equipment supplied against the delivery note. Please report incomplete deliveries immediately to the shipping carrier or the supplier.

Interim Storage

The external packaging of the unit and its spare and replacement parts is designed for a storage period of three months from delivery.

Units kept in interim storage must be stored in a closed, dry room at temperatures from +5°C to max. +50°C.

5.3 Installation Instructions

5.3.1 Ambient Conditions and Conditions at the Place of Installation

Select a place of installation that meets the following conditions:

- ▶ No relatively high temperature or exposure to heat or sun
- ▶ No strong drafts
- ▶ No aggressive chemical environments
- ▶ No strong vibrations
- ▶ No dust
- ▶ No major humidity effects

5.3.2 Transporting and Unpacking the Equipment

The unit and its accessories are packaged according to the national custom and protected against impacts. To avoid damaging the equipment, please comply with the following notes and instructions:

- ▶ Please observe the Technical Specifications in Chapter 3!
- ▶ Please observe the weight. For transport to the unpacking site, use a floor truck if appropriate.
- ▶ To prevent the unit from slipping, it should always be carried or lifted by two people.
- ▶ To open the packaging, remove the unit and place it on a workbench or table intended for that purpose.
- ▶ The workbench or table must be level and able to bear the appropriate load. Allow sufficient space around the workbench or table for the operator to work.
- ▶ Unpack the accessories and technical documents and put them in a safe place.
- ▶ Moisture in the air can condense on the surfaces of a cold unit whenever it is brought into a substantially warmer place. If you transfer the unit to a warmer area, make sure to condition it for about 2 hours at room temperature, leaving it plugged into AC power. The positive temperature differential between the inside of the unit and the ambient environment essentially protects the unit against humidity and moisture.

Open the packaging carefully. It is a good idea to save the packaging material (in case the carrier wants to inspect it or the shipment has to be returned).

Dispose of packaging materials in compliance with prevailing environmental safety regulations.

5.4 Connections

5.4.1 Front Panel and Back Panel

Fig. 5-2 | Front Panel


Fig. 5-3 | Back Panel

5.4.2 Left and Right Panels Sartocheck® 4 plus Filter tester


Fig. 5-4 | Left panel

Before starting up the unit; make sure that the unit is stable and secure. Check all pneumatic and electrical connections on the unit.

The following displays are located on the front panel:

- 1 Standby**
Is lit up when Sartocheck® 4 plus Filter tester is ready for operations, but not currently running any test program. (Blinking LED)
- 2 Ready**
Lights up when the result of the filter test is ready. (LED lit)
- 3 Running**
Is lit up when a test is currently running. (LED lit)
- 4 Error**
Lights up when the test unit has identified an error and cancelled the test. (LED lit)

The following connection ports are located on the back panel:

- 1 Housing ventilator**
The unit is equipped with a ventilator that ensures that Sartocheck® 4 plus Filter tester is kept at a constant temperature.
- 2 Serial port MU (Management Unit)**
This port is used for connecting accessories (e.g. barcode scanner).
- 3 AC adapter**
Power plug for connecting the power cord to mains supply.
- 4 Fuses**
Fine-wire fuses between unit and mains supply.
- 5 Network Connection**
RJ45

Located on the left panel:

- 1 Main power switch**
Switch on the power with the main power switch.
- 2 SD Card Reader**
The SD card reader is designed for standard SD cards for inputting and storing data.
- 3 Serial port TU (test unit)**
Port for connecting additional test units (Multiunits).
- 4 PLC Port**
Port for connecting PLC controller signal input and output.


Fig. 5-5 | Right panel

Located on the right panel:

- 1 Ext. sensors**
Port for connecting an external pressure sensor.
- 2 Ext. valves**
Port for connecting external valves.
- 3 Ext. reference tank**
Quick connector for external reference tank.
- 4 Venting 1**
Quick connector for Venting 1.
- 5 OUT**
Quick connector for test tubing.
- 6 Venting 2**
Quick connector for Venting 2.
- 7 IN**
Quick connector for compressed air supply.

5.4.3 Left and Right Panels Sartocheck® 4 Filter tester MultiUnit

Fig. 5-6 | Left panel

Fig. 5-7 | Right panel

5.5 Pneumatic Connections

Sartocheck® 4 plus can be operated with an internal pressure sensor and with different ports for connecting tubing. The ports for connecting tubing, for both the Sartocheck® 4 plus and Sartocheck® 4 MultiUnit, are illustrated in Fig. 5-8, Fig. 5-9 and Fig. 5-10.

Do not connect Sartocheck® 4 plus Bag tester to pressurized or high temperature systems! Wait until the whole system is at room pressure and temperature!

Inlet Pressure During Operation with Internal Pressure Sensor

The minimum operating pressure required depends on the test method selected, but should be at least:

> 4000 mbar and > 500 mbar above the test pressure or above the max. B.P. entered. Sartocheck® 4 plus automatically detects the operating pressure during the function test before starting each test and during the whole test sequence.

Check to make sure that all inlet and outlet pressure tubing is connected properly, that all electrical connections have been correctly installed, and that the paper roll as well as the ribbon cartridge have been inserted.


Fig. 5-8 | Connection with internal pressure sensor


Fig. 5-9 | Connection with internal pressure sensor and external vent


Fig. 5-10 | Connection with internal pressure sensor and external volume reference tank

5.6 Connections with External Pressure Sensor

Sartocheck® 4 plus can be operated with an external pressure sensor and with different ports for connecting tubing. The ports for connecting tubing, for both the Sartocheck® 4 plus and Sartocheck® 4 MultiUnit, are illustrated in Fig. 5-11, Fig. 5-12 and Fig. 5-13.

Check to make sure that all inlet and outlet pressure tubing is connected properly, that all electrical connections have been correctly installed, and that the paper roll as well as the ribbon cartridge have been inserted.


Fig. 5-11 | Connection with external pressure sensor and manual valves


Fig. 5-12 | Connection with external pressure sensor and external valves


Fig. 5-13 | Connection with external pressure sensor for WIT

5.6.1 Inserting Paper Roll Ribbon Cartridge

To insert or change the paper roll, lift the cover off the paper feed. Pull the feed mechanism out of the assembly and push the cover to the left and lift it out of its attachments. Now lift the paper feed mechanism out of the Sartocheck® 4 plus housing.

Insert a new paper roll into the paper slot and slide the paper feed mechanism back into the housing. Thread the end of the paper through the slot of the paper feed.

Make a final check to ensure that the feed mechanism is locked in place.

Close the printer cover and secure the splashguard. Check to make sure that the splashguard is locked in place.

To change the ribbon cartridge (when the feed mechanism is pulled outwards) lift the previously installed cartridge out of the paper feed and insert a new one. Apply slight pressure to lock the cartridge in.

5.6.2 Sartocheck® 4 plus
Filter tester MultiUnits

Sartocheck® 4 Bag tester MultiUnits are connected to the Sartocheck® 4 plus Filter tester MultiUnit or to the Sartocheck® 4 plus Bag tester via the serial port TU connection ports. Different ways of connecting are shown below in Fig. 5-14 to 5-17. The location of the Sartocheck® 4 Bag tester MultiUnits, relative to the Sartocheck® 4 plus Filter tester or to the Sartocheck® 4 plus Filter tester MultiUnit, will determine which method is best suited for the application.


Fig. 5-14 | Serial connection to a Bag tester


Fig. 5-15 | Serial connection to a Filter tester

6. Function

Fig. 6-1 | Overview

Picture of Sartocheck® 4 plus Filter tester

- 1 Housing
- 2 Printer
- 3 Touch screen
- 4 Keypad
- 5 Indicator lights

6.1 General Functional Description

Sartocheck® 4 plus Filter tester is a qualifiable filter integrity testing unit, consisting of a data management unit, a test unit and up to four additional Filter tester MultiUnits or Bag tester MultiUnits.

The test unit is a mechanically, pneumatically and electrically closed functional unit.

A test unit cannot be operated directly by the user, but is 100% controlled by the data management unit via an internal RS 485 interface. With the help of its electrical and pneumatic peripheral devices, the data management unit performs the actual filter tests along with a series of other additional functions.

Bag Test

The main purpose of Sartocheck® 4 plus Bag tester is to test single use bags and single use bioreactors according to the pressure drop method.

Filter Integrity Test

The main purpose of Sartocheck® 4 plus Filter tester is to test the integrity of membrane filters according to various methods.

6.2 Test Programs

In the test program, the procedural steps and the parameters are set down for either the bag test or for the filter integrity test.

A programmed test program can be loaded into the test unit and the test will be run independently.

- Load a test program
- Run a test program

The test unit runs through the previously loaded test program, in other words, it carries out the test and saves the result in non-volatile memory. The data management unit monitors the progress of the bag test or the filter integrity test, i.e. it reads in the current data in order to present them to the user.

6.3 Managing the Test Results	Sartocheck® 4 plus Filter tester processes the measured raw data (the test results), which means that the unit values are converted and numbers transformed into text.
6.4 Function Test	With the function test, the user is able to interactively test the functionality of pneumatic components. The proportional valves are not explicitly, but rather implicitly tested in this context.
6.5 Cleaning	Sartocheck® 4 plus Filter tester starts the procedure for cleaning and flushing the internal lines and the volume reference tank.
6.6 User Interrupt	The user can interrupt any test on Sartocheck® 4 plus Filter tester by pressing the [STOP] button. Pressing the [arrow upwards] and [STOP] buttons simultaneously will interrupt tests on all Sartocheck® 4 plus Filter tester MultiUnits or Sartocheck® 4 Bag tester MultiUnits connected to the Sartocheck® 4 plus Filter tester, as well as any test running on the Sartocheck® 4 plus Filter tester unit.
6.7 Self-Test for Broken Conductors and Short-Circuiting	While controlling the valves and sensors, the Sartocheck® 4 plus Filter tester hardware can check whether there are any broken conductors or short-circuits.
6.8 Calibration	The calibration function is designed to calibrate the internal pressure sensor with the help of more accurate reference measuring equipment, and if necessary, adjust it. The sensor has a monitored calibration date. An adjustment routine automatically follows the test. The test can also be activated separately. Sartocheck® 4 plus Filter tester saves the raw data for the calibration record in permanent memory. The operator can use the following sensor type for determining the reference pressure: ► Relative pressure sensor: The measured reference pressures relate to the ambient pressure. The ambient pressure must be known or measured with additional reference measuring equipment.
6.9 Internal Pressure Sensor	Sartocheck® 4 plus Filter tester sets 8 pressure values and compares them with the internal pressure sensor and the pressure values measured on the reference measuring instrument (connected to the "OUT" connector) to calculate scale and linearization.
6.10 Pressure Gauge Function	When in the pressure gauge function, Sartocheck® 4 plus Filter tester acts like a digital pressure measuring device.
6.11 Cleaning Drying	Through the pneumatic connectors, Sartocheck® 4 plus Filter tester is able to flush almost all internal volumes. It is the responsibility of the user to use original Sartocheck® equipment only (Cleaning Kit type 26288---CK).
6.12 Safety Functions	► As one of its technical options, Sartocheck® 4 plus Filter tester runs the test step Initialization Security Status in the event of a power failure. Basically, what this step does is close the compressed air supply and the block valve to the test set-up. ► Sartocheck® 4 plus Filter tester is not responsible for ensuring whether the pressure supplied to the system being tested is too high. This means that the user must make sure that the programmed pressure cannot damage the system being tested.

6.13 Software Update

The operating software on Sartocheck® 4 plus Filter tester is updated with a corresponding program run with an external PC | laptop. The Management Unit may be updated with an SD card. Since the configuration data are localized in a battery-backed RAM, they remain unchanged during any software update.

6.14 Serial Number

Every unit bears its own unique serial number. The serial number may not be changed during a software update.

- ▶ The serial number is entered for the first time during production of the unit.
- ▶ If the PCB has to be replaced during repairs that may arise later, the unit must be assigned the exact same serial number that is imprinted on the manufacturer's label.

7. Operation

7.1 Switching On the Unit

7. Operation

Fig. 7-1 | Overview

Picture of Sartocheck® 4 plus Filter tester

- 1 Housing
- 2 Printer
- 3 Touch screen
- 4 Keypad
- 5 Indicator lights

7.1.1 Operating Elements


Fig. 7-2 | Display and control panel

The display and control panel features the following operating elements:

- 1 Esc**
Press this key to cancel the current procedure.
- 2 F1-F 10**
Function keys F1-F10 correspond to the function types displayed on touch screen.
- 3 Log on**
Press this key to switch on the log-on screen.
- 4 Log off**
Press this key to log the current user off.

5 Start

Press this key to start a test.

6 Stop

Press this key to stop a running test.

7 Paper

This key advances the paper in the internal printer.

8 PC keypad

Is used for inputting and starting processes.

9 Home

When this key is pressed, programming is cancelled without saving. The former test is reinstated and the user is automatically taken back to the Main Menu.

10 End

Press this key to jump to the end of the input line when entering parameters and data.

11 Page up

Press this key to scroll a page back in the service menu guidance.

12 Page down (toggle key)

Press this key to scroll to the next service menu window.

13 All Test Stop

Pressing the [↑] and [STOP] buttons simultaneously will stop all running tests on all MultiUnits, as well as any running test on the Sartocheck® 4 plus Filter tester unit.

Both the keypad and the touch screen can be used to operate Switching Sartocheck® 4 plus Filter tester. These two options can also be combined.

Individual menu items can be selected Sartocheck® 4 plus Filter tester with function keys, on the keypad or by pressing the touch screen. Selection from lists is done by entering the line number or by direct selection on the touch screen.

All dialog boxes are set up so that the operator can tell by the color which part of the touch screen can be operated, which are the input boxes and which areas are for information only. Function keys and other keys that can be pressed on the screen and on the keypad are highlighted in yellow.

Input boxes are white. Non-changeable information boxes or function keys that have already been activated are highlighted in gray. Selected dialog boxes are colored blue. Warnings and error messages appear in red.

7.2 Switching the Unit
ON and OFF

If you want to switch the unit on, proceed as follows:

Before you switch on the power supply and the compressed air, compare the connection data for the unit (see the Chapter 01 Technical Specifications) with the data for the mains current and the compressed air supply.

Connect Sartocheck® 4 plus Filter tester to the mains supply by plugging the power cord into the wall outlet, install all tubing connections. The unit is now ready for operation.

Fig. 7-3 | Login menu

Login menu

User ID:	Enter operator's name
Password:	Enter operator's password
OK:	Confirm correct entry
Cancel:	Delete entry
Important Note:	When Sartocheck® 4 plus Filter tester is delivered it is factory set for the following user:
Name:	mentor
Password:	mentor

Operators must log-in before activating the Main Menu. The operator logs in by pressing the "Login" key on the keypad. After pressing the "Login" key, the user will be asked for their name and password. When both are entered correctly, the keypad and touch screen are released for use. When the operator stops working with Sartocheck® 4 plus Filter tester, they must always log out by pressing the "Logout" key. Sartocheck® 4 plus Filter tester is then blocked for further operation until another operator logs in.


Fig. 7-4 | Main Menu

Information line with model, version, system number and user name.

Menu display: Main Menu

F1-F8: Select the function menus with the function keys

Selects the function menus

Functional unit selected and operational status.

The Main Menu appears.

- Before starting a test, allow the unit to warm up for about 20–30 min. For shut-down times < 30 min and constant ambient temperatures, the warm-up phase can be shortened accordingly. For more information about the tests, see Chapter 7.3 ff.
- Program the data for the test routine.
- The Sartocheck® 4 plus Filter tester ready for operation (by the logged-in user).
- Start the test.

Only service technicians or trained personal may install connections and start up the unit!

7.2.1 Running a Test

Run the test according to the test procedure that you selected and programmed.

7.2.2 Switching Off the Unit

When switching off the unit, proceed as follows:

- ▶ Run the test to the end.
- ▶ Close the compressed air supply at the compressed air supply port.
- ▶ Run the log-off function.
- ▶ Turn off the main power switch.

7.2.3 Switching Off the Unit in the Event of a Power Failure

In the event of a power failure, all valves will be switched to a secure state. This guarantees that the system being tested is vented. The ongoing test is cancelled. No test result will be available.

Loss of Compressed Air

In the event that the compressed air supply stops, the test is automatically cancelled and an error message appears. "Target pressure not reached". The system is vented automatically. When the supply of compressed air resumes, the unit must be restarted and the test repeated.

7.3 Tests

7.3.1 Test Parameters

To obtain reliable and reproducible test results, observe the following:

- ▶ Run the test at room temperature (approx. 15–30°C), whenever possible.
- ▶ During a test run, always avoid any temperature changes.
- ▶ The upstream side of the system must be absolutely leak-free
- ▶ Do not connect the tester to pressurized systems.
- ▶ For filter testing: If you wish to use other types of compressed gas and/or wetting fluids not specified by Sartorius Stedim Biotech, consult your local Sartorius Stedim Biotech office or Sartorius Stedim Biotech GmbH in Germany.)
- ▶ For the Water Intrusion Test, the surface tension of the water must be > 70 dyn/cm.
- ▶ Always enter correct test parameters. Incorrect test parameters will lead to incorrect test results or erroneous evaluations.

Calibration

The Sartocheck 4 plus Filter tester should be re-calibrated at least once a year, twice being preferable.

7.3.2 Test Methods

With the Sartochek[®] 4 plus Filter tester, you can directly select a particular test method. This greatly simplifies programming the data for a test, since only the parameters that are relevant for the selected test are displayed.

You can select any of the following test methods:

Test methods

- 1 Diffusion test
- 2 Bubble point test (Standard and Customer specific)
- 3 Diffusion and bubble point test
- 4 Water intrusion test
- 5 Pressure drop test
- 6 Multipoint diffusion test (Standard and Customer specific)
- 7 Volume measurement
- 8 Plate holder test
- 9 Installed bag test

System Venting

At the end of each test, after a power failure or after switching off the unit during a test, the system is automatically vented.

7.3.3 Programming Tests
(F2 – Main Menu)

Fig. 7-5 | Main Menu

The Sartochek[®] 4 plus Filter tester cannot perform bag tests. Nevertheless it is possible to program bag tests due to the fact that additional Bag tester MultiUnits may be connected to the Filer tester. The test parameters for the individual tests are entered and selected on Sartochek[®] 4 plus Filter tester by using the keypad or the touch screen. The Main Menu is the starting point for programming. The following example uses a diffusion test to illustrate how to program data for tests and run test sequences (F2 – Main Menu). The programming and test sequences are identical for the other test types.

Before we start programming, we would like to explain a few important keys and/or touch commands.

The Sartochek[®] displays are set up with yellow buttons located on the right side of the screen. These buttons always have the same functions:

- Main Menu: The program jumps back to the Main Menu (see Fig. 7-5)
- Page Up: the program scrolls to the previous page of the screen
- Page Down: the program scrolls to the next page of the screen
- List Up: Whenever a list is displayed on the screen, use LIST UP to scroll up the list
- List Down: Whenever a list is displayed on the screen, use LIST DOWN to scroll down the list
- Enter: Confirms the selected function, action or entry. Following confirmation, the yellow key turns gray


Fig. 7-6 | Programming tests: Program selection

Programming Tests

Select the programs in the menu "Programming tests: Program selection". You can select existing programs by entering the line numbers or by selecting them directly on the touch screen. The selected line is highlighted in blue. Scroll back and forth in the program selection list with List DOWN | UP. Confirm the selected program line by pressing Enter. The Enter key turns gray. If you would like to correct the program memory location, select a new line number or highlight the new line directly on the touch screen.

Re-confirm by pressing Enter. Press PAGE DOWN to scroll to the next page. The diffusion test is used as an example.

If the optional Barcode Scanner is used for programming, please compare chapter 7.18 on page 49).


Fig. 7-6b | Optional parameters

Optional Parameters for Detection of Disconnected Filters

In case of a disconnected filter (housing), an integrity test can be successfully performed because the test-evaluating parameters (e.g. BP_{min} , $Diff_{max}$) are fulfilled. To avoid such situations, it is possible to define four additional parameters. A non-standard test (compare fig. 7-8) must be chosen to get access to these parameters. As soon as at least one of these listed parameters are not fulfilled, a warning message is displayed | printed. The following parameters can be defined:

Min. limits: Parameters can be activated | deactivated
The parameters can be activated / deactivated.
If deactivated the here below described limits will have no impact on the test sequence.

Min. DIF/WIT/WFT: A minimal diffusion rate can be defined.
This parameter is valid for diffusion test, WIT and WFT.
If the measured value falls below this limit a warning message will be printed at the end of the test together with the test evaluation. It is then responsibility of the of the quality assurance to decide if the test is conforming or not.
The test evaluation of the Sartocheck 4 plus remains "passed".

Min. Flow BPT: A minimum flow rate can be defined for the Bubble Point Test.
This value represents the minimum flow that must be measured at BP_{max} .
If the measured value falls below this limit a warning message will be printed at the end of the test together with the test evaluation. It is then responsibility of the quality assurance to decide if the test is conforming or not.
The test evaluation of the Sartocheck 4 plus remains "passed".

Min. Volume: A minimal net volume can be defined. A net volume measurement below this value will abort the test. An error message will be printed.

Max. Volume: A maximal net volume can be defined. A net volume measurement above this value will abort the test. An error message will be printed.

7.4 Diffusion Test

Fig. 7-7 | Programming tests: Program selection

Fig. 7-8 | Test parameters: External connections

When pressure is applied to the upstream side of a wet filter, a diffuse gas stream starts to flow through the filter membrane after a short time. Diffusion is dependent on the effective filter surface area. The diffusion test is not suitable for testing disk filter systems. However, the diffusion test is a very well suited test method for the integrity testing of medium-sized and larger filter cartridge systems starting from approx. 150 cm². The diffusion value refers to a test pressure, which means there is a direct correlation to the actual retention behavior of the filter (for a more detailed description, see Chapter 12, Appendix).

In the menu "Programming tests: Program selection", the test type is selected.

- Select test type (standard test or user-defined test)

In the menu "Test parameters: External connections," decisions about external connections and the configuration are entered. The function keys are used here to enter the appropriate settings. The following external connections and parameters must be configured for all tests:

- Standard test (Yes|No)
- External pressure sensor (Yes|No)
- External valves (Yes|No)

The parameters for the standard test are defined in the configured Service Menu. Unlike the defaults (standard values), the values can be individually adjusted.

For the diffusion test, you have the following options:

Pressure unit: mbar | psi | hPa
Automatic test mode: No | Yes

Enter parameters for the diffusion test

Test pressure: from 100 – 8000 mbar
Stabilization time: in min.
Test time: in min.
Max. diff.: in ml/min
Net volume: in ml


Fig. 7-9 | Test parameters: Diffusion test

In the menu "Test parameters: Diffusion test," the parameters for the diffusion test are entered with or without net volume measurement.

For tests with volume measurements:

- Test pressure
- Stabilization time
- Test time
- Max. diffusion

When conducting tests without volume measurement, net volume should be entered separately.

The test parameter data can be found in the validation documents of the respective filter cartridge.

Locks

The locks are activated for the test parameters (default). After entering with the Shift + F keys, the programmer can lock the data to prevent any changes. The dialog boxes are then highlighted in gray. The user can no longer change the entered data.

Fig. 7-10 | Test program: Record data 1

Enter record data 1 (up to 32 characters per field)
 Company: Company name
 Building: Building name; house number
 Department: Identification of the department
 Manufact. site: Production sector
 Product: Product name
 Product lot: Batch number
 Filter: Filter size
 Filter lot: Number of the filter batch

Fig. 7-11 | Test program: Record data 2

Locks

After entering with the Shift + F keys, the programmer can lock the data to prevent any changes. The dialog boxes are then highlighted in gray. The user can no longer change the entered data. The relevant record data are entered in the menu "Test program: Record data 1."

Entry of record data 2

Filter line: Name of filter line
 Housing: Type of housing
 Wetting medium: Type of wetting medium
 Test gas: Type of test gas
 Water quality: Data on water quality
 Comment: Name of operational process
 Comment 2: Enter test-relevant remarks
 Comment 3: Enter test-relevant remarks



Fig. 7-12 | Test program: Start

Start the Test Program

Save:
 After confirmation, the data are saved

The menu "Test program: Start" lists the options for that test start.

Once the programming is finished, the test can be started.


Fig. 7-13 | Test unit selection

Choose the test unit on which to perform the integrity test.
 Only test units that are currently available will be selectable.


Fig. 7-14 | Test: Start

The test sequence is structured as follows:

Pressurization

Sartochek® 4 plus Filter tester adjusts the pressure to the pre-programmed test pressure.

After a short waiting period during which the system is tested for major leakage, the program goes into the stabilization phase.

Fig. 7-15 | Stabilization

Stabilization Phase

During the stabilization phase, the test pressure is adjusted to the exact desired value. Fluctuations in test pressure occur as a result of cartridge compaction and temperature variations created by the different filter housing temperatures and the downstream air flowing in from the compressed air supply.

These fluctuations in test pressure are compensated for during stabilization. About half of the stabilization time is utilized to carry out volume measurement (when volume measurement is selected).

Test Phase

In the test phase, the system is in a closed state. The pressure drop resulting from diffusion is measured over the test time entered.

Fig. 7-16 | Volume measurement

Volume Measurement

A clock is shown during the volume measurement.


Fig. 7-17 | Volume Measurement

Volume Measurement
The volume is displayed in ml.


Fig. 7-18 | Test Phase

In the test phase, the current diffusion value is displayed as a graph and additionally as a numerical value. The graph is designed so that the limit is drawn as a red line. If the actual value is below this, the integrity of the filter is proven.

7.5 Bubble Point Test

In wet membranes, the wetting liquid is held inside the porous structure by capillary force. This force increases as the pore size decreases. A specific gas pressure dependent on the pore size of the membrane is necessary to force out the liquid from the pores. This pressure is generally designated as the "Bubble Point," or B.P. for short. It generally indicates the largest pores of a filter since liquid is first expelled from them. Hence, it follows that the B.P. of a filter depends both on the filter material (the contact angle varies from material to material) and, in particular, on the surface tension of the wetting fluid. Low surface tensions, as are present in organic solvents, detergents and emulsifiers yield low B.P.'s.

Since the surface tension also depends on the temperature, different wetting fluid temperatures will also yield different B.P.'s. The Bubble Point Test is especially suitable for checking the integrity of disc filter systems. (For a more detailed description, see Chapter 12, Appendix).

Test parameters: Bubble Point Test


Fig. 7-19 | Test parameters: Bubble Point Test

Enter the parameters: BP Min. in mbar
BP Max. in mbar
System test class: e.g. standard
Net volume: Enter or measure, if Enter: Volume in ml

Recommendation for test-class selection:
Small Systems: disc filters, mini cartridges, capsules up to 5"
Standard Systems: 10"–20" filter cartridges and capsules
Special Systems: 30" and small multiple cartridge systems

For bigger systems please use the customer specific BP Test.
Please contact Sartorius Stedim Biotech if specific settings are needed.

7.6 Customer-Specific BP

The normal bubble point test (as previously described) uses a standard set-up of parameters to cover most of applications. However, some applications in which the wetting fluid differs considerably from water (e.g. oily or viscous products) or in which the filtration area is large (multi-round housings), bubble point detection can be optimized using the customer-specific bubble point test. Please ask for Sartorius Stedim Biotech support to program the optimal parameters.

Test parameters: customer-specific BP

- Start factor:
- Enter the factor
- Prestabilization:
- in sec.
- Stability by phase:
- in sec.
- Phase test:
- in sec.
- A1 Criterion :
- in ml/min
- A2 Criterion :
- in ml/min


Fig. 7-20 | Test parameters: Customer-specific BP

7.7 Diffusion and Bubble Point Test (Combined Test)

This test method is a combination of the diffusion and the B.P. tests and is predominantly used to check the integrity of membrane filters.

Test parameters: Diffusion and Bubble Point Test

- Test pressure:
- in mbar
- Stabilization time:
- in min.
- Test time:
- in min.
- Max. diffusion:
- in ml/min.
- Min BP:
- in mbar
- Max BP:
- in mbar
- Net volume:
- Enter or measure, if Enter: Volume in ml


Fig. 7-21 | Test parameters: Diffusion + BP Test

7.8 Water Intrusion Test |
Water Flow Test

This test method is used exclusively to test hydrophobic gas filter elements. The test sequence of this test method is the same as that of the Diffusion Test. However, the Water Intrusion Test differs from the latter test in that the hydrophobic filter element in the filter housing is flooded upstream with water to blind the filter so that when a pressure gradient is applied, the water volume intruding into and diffusing through the hydrophobic filter matrix per unit of time is calculated from the pressure drop, instead of the gas stream diffusing through the membrane as in the Diffusion Test.

- 1. The Water Intrusion Rate: Indicates the rate calculated from the pressure drop, test time and upstream volume (analogous to the calculation of diffusion).
- 2. The Water Flow Rate: Indicates the rate related to the water flow in the membrane (Water Intrusion converted to "Water Flow").

The Water Intrusion Test is suitable for all hydrophobic systems of approx. 1000 cm² filtration area and larger, all the way to process-scale systems. (For a more detailed description, see Chapter 12, Appendix)

Test parameters: Water Intrusion Test

Field	Value	Unit
Program no.	4	
Name	WITdemo	
F1 Test pressure	2500	mbar
F2 Stabil. time 1	3	min
F3 Stabil. time 2	10	min
F4 Test time	10	min
F5 Max. WIT	13.0	ml/10 min
F6 Auto. filling	No	
F7 Net volume	Measurin	
F8 Volume	0	ml

Fig. 7-22 | Test parameters: Water Intrusion Test

- Test pressure: in mbar
- Stabil. time 1: in min.
- Stabil. time 2: in min.
- Test time: in min.
- Max. WIT: in ml/min
- Auto filling: Yes| No
- Net volume: Enter or measure, if Enter: Volume in ml

Remark:
Stab 1 is the stabilization time before Netvolume determination. This time is needed to compact the filter system and therefore to avoid misinterpretation. Stab 1 should not be less than 3 min.
Stab 2 is the stabilization time before the water-intrusion (water-flow) measurement. The stab 2 value is available in the validation guide of the filter cartridges. The standard value is 10 min.

7.9 Pressure Drop Test

The Pressure Drop Test is based on the change in pressure within a filter system connected upstream. This change is caused by the diffusive gas stream flowing through a wet membrane at a given test pressure (see "Diffusion Test"). The pressure drop measured in the process can be used as a direct measure of the filter system's integrity, provided that the upstream volume of the system is constant and known. For this reason, a Pressure Drop Test is recommended mainly to leak-test empty housings and pressure vessels with a known volume.

We additionally recommend that you carry out the following Pressure Drop Test to check for major leakages. (For a more detailed description, see Chapter 12, Appendix).

Pressure Drop Test

Test pressure: 3000 mbar
Stabilization time: 3 min
Test time: 3 min
Pressure drop: max. 3 mbar

Test parameters: Pressure Drop

Test pressure: in mbar
Stabilization time: in min.
Test time: in min.
Max. pressure drop: in mbar


Fig. 7-23 | Test parameters: Pressure Drop Test

7.10 Multipoint Diffusion Test

The Multipoint Diffusion Test is based on the dependence of the diffusive gas flow to the actual test pressure. Starting at a given minimum pressure, the test pressure is increased step by step up to the maximum value, or the bubble point, and the diffusion is measured for each pressure increment. All values are output by Sartocheck® 4 plus in the form of a pressure | diffusion curve and a data table.

The slope of the curve allows changes in the filter material to be identified before the next filtration run so that by comparing the curves obtained at the beginning of the filter's use, any defects can be detected before they can spoil an actual filtration run. The Multipoint Diffusion Test is thus suitable as an additional test method for verifying integrity test results.

Test parameters: Multipoint Diffusion Test

Min. test press.: in mbar
Max. test press.: in mbar
Pressure incr.: e.g. in 100 mbar increments
Test class: e.g. Standard
Net volume: Enter or measure, if Enter: Volume in ml


Fig. 7-24 | Test parameters: Multipoint Diffusion Test

7.11 Customer-Specific
Multipoint Diffusion Test

Support point	mbar	ml/min
Support point 1	1100	10.0
Support point 2	1200	20.0
Support point 3	0	0.0
Support point 4	0	0.0
Support point 5	0	0.0
Support point 6	0	0.0
Support point 7	0	0.0
Support point 8	0	0.0

Fig. 7-25 | Limit line parameters

The normal Multipoint Diffusion Test does not determine a pass or fail result. This Customer-Specific Multipoint Diffusion Test allows the user to enter a maximum diffusion limit line. Up to eight (8) points can be used to define this limit line (pressure vs. diffusion). During test execution the actual measured diffusion will be compared to this limit line to determine a pass or fail result.

7.12 Volume Measurement

This test method measures the net volume of the system (e.g. filter housing). The Sartochek® 4 plus Filter tester uses its internal or external reference tank pressurized to the programmed pressure to determine the unknown volume. This is not an integrity test. Rather, this test enables net volume validation of vent filter installation for water intrusion tests using an external pressure sensor, for instance. Please ask for Sartorius Stedim Biotech support before performing any net volume validation.

Enter test parameters: Volume Measurement

Test parameter	Value
Test pressure	1000 mbar
Stabilization time	1 min

Fig. 7-26 | Test parameters: Volume Measurement

Test pressure: in mbar
Stabilization time: in min.

7.13 Plate Holder Test

This test sequence may be programmed on the Sartochek 4 plus Filter tester but may only be performed if a Bag tester MultiUnit is connected. The test sequence is designed for leak testing of 2 dimensional bags (e.g. Flexboy) being maintained between 2 metal plates during the test e.g. using the FlexAct BT trolley equipped with fleeces. The test sequence is based on the change in pressure within the bag system including connected tubes.

The pressure drop measured in the process can be used as a direct measure of the bag set-up leak tightness and absence of pinholes or similar defects, provided that the upstream volume of the system is constant.

Enter parameters for the plate holder test

- Test pressure: from 3–300 mbar
- Stabilization time: 10-3600 seconds.
- Test time: 10-3600 seconds
- Max. pressure drop: 1-programmed test pressure

Please verify the test parameters for the bag being tested. Applying too high a test pressure may jeopardize the integrity of the bag.

Since the viscosity of the remaining liquid may influence the leak rate of an eventual pinhole, please do not mix up pre-use test values with specifically validated post-use values.

This test can also be used for testing vessels as long as they resist the applied pressure.

Enter record data 1 (up to 32 characters per field)

- Company: Company name
- Building: Building name; house number
- Department: Identification of the department
- Manufact. site: Production sector
- Product: Product name
- Product lot: Batch number
- Bag type: E.g. Flexboy 50L
- Bag lot: Number of the bag batch


Fig. 7-27 | Test parameters for plate holder test

Fig. 7-28 | Data log group 1

Locks

After entering with the Shift + F keys, the programmer can lock the data to prevent any changes. The dialog boxes are then highlighted in gray. The user can no longer change the entered data. The relevant record data are entered in the menu "Test program: Record data 1."

Entry of record data 2

- Nom. bag volume: Nominal volume of the bag
- Test gas: Type of gas
- Comment: Name of operational process
- Comment 2: Enter test-relevant remarks
- Comment 3: Enter test-relevant remarks

Fig. 7-29 | Data log group 2

7.14 Installed bag Test

This test sequence may be programmed on the Sartocheck 4 plus Filter tester but may only be performed if a Bag tester MultiUnit is connected. The test sequence is designed for leak testing of 3 dimensional bags (e.g. STR Cultibags) being installed in the appropriate STR holder equipped with fleeces.

The test sequence is based on the change in pressure within the bag system including connected tubes and vent filters.

The pressure drop measured in the process can be used as a direct measure of the bag set-up leak tightness and absence of flaws or similar defects, provided that the upstream volume of the system is constant.

Sartocheck 4plus V.YXX S/N:12345678 User name: Mentor

Test parameters Installed bag test Main Menu

Program no.: 2 Name: Joli

F1	Test pressure:	30.0 mbar	3	Page up
F2	Stabilization time:	1 min	3	
F3	Test time:	2 min	3	Page down
F4	Max. pressure drop:	5.0 mbar	3	
F5	Filling time:	5 min	3	

Enter

Bag test unit 1: MultiBag Test unit 2: OffTest unit 3: OffTest unit 4: OffTest unit 5: Off

Fig. 7-30 | Test parameters for installed bag test

Enter parameters for the plate holder test

Test pressure: from 3–50 mbar
Stabilization time: 1–1080 minutes.
Test time: 1–1080 minutes
Max. pressure drop: 1-programmed test pressure
Filling time: 1–1080 minutes

Please verify the test parameters for the bag being tested. Applying too high a test pressure may jeopardize the integrity of the bag.

Since the viscosity of remaining liquid may influence the leak rate of an eventual flaw, please do not mix up pre-use test values with specifically validated post-use values.

This test can also be used for testing vessels as long as they resist the applied pressure.

Enter record data 1 (up to 32 characters per field)

Company: Company name
Building: Building name; house number
Department: Identification of the department
Manufact. site: Production sector
Product: Product name
Product lot: Batch number
Bag type: E.g. STR Cultibag 1000L
Bag lot: Number of the bag batch

Locks
After entering with the Shift + F keys, the programmer can lock the data to prevent any changes. The dialog boxes are then highlighted in gray. The user can no longer change the entered data. The relevant record data are entered in the menu "Test program: Record data 1."

Entry of record data 2

Nom. bag volume: Nominal volume of the bag
Test gas: Type of gas
Comment: Name of operational process
Comment 2: Enter test-relevant remarks
Comment 3: Enter test-relevant remarks

7.15 Program Database (F4 – Main Menu)


Fig. 7-31 | Program database

All test programs created are saved in the program database.

Sartochek® 4 plus Filter tester features the option to archive 250 user-programmable test programs in its internal database. The program data in the program item "Program database (F4)" can be managed.

The following functions can be selected:

- F1: Load test programs from internal memory
This will display all programs inside the Sartochek® 4 plus Filter tester
- F2: Upload test programs from SD card
This gives the possibility to install new programs in the Sartochek® 4 plus Filter tester (e.g. transfer from one Sartochek® 4 plus to another)
- F3: Upload test program from network
Allows uploading test programs from the network, eventually previously downloaded from another Sartochek® 4 plus, if the if FTP directory is identical for both Sartochek's.
- F4: Save backup to SD card
Saves a full copy of the programs set-up to an SD card
- F5: Save backup to network
Saves a full copy of the programs set-up to the network
- F6: Load backup from SD card
Restores a copy of test programs set-up from an SD card (this will overwrite the existing programs)
- F7: Load backup from network
Restores a copy of test programs set-up from the network (this will overwrite the existing programs)
- F8: Delete PLC programs

With F4 and F5, the test programs are stored in backup format.

This applies to the entire program database only. Individual programs cannot be stored in backup format.

7.15.1 Loading Test Program from Internal Memory (F1)


Fig. 7-32 | Program database: Internal memory

When loading test programs from internal memory, all test programs are displayed in the program selection list. You can highlight the programs by entering the list number or using the touch screen. There is an entry line "From ## to ##" above the program list. In this line, you can highlight test program blocks.

Confirm your selection by pressing enter.

The highlighting can be removed with the button "Undo marking".


Fig. 7-33 | Program database

The highlighted test program is displayed above the function keys F1-F4.

When test program blocks are highlighted, you can scroll through the selection list with LIST UP|DOWN.

You can select the following functions:

- F1: View
Display the parameters predefined in the test program along with previously defined record data.
- F2: Print
Print out the test program on the internal printer
- F3: Save parameters to SD card
Save the test program to an SD card
- F4: Save to network
Save the test program to the FTP server
- F5: Save as program PLC 1
Save the selected program to memory location PLC 1
- F6: Save as program PLC 2
Save the selected program to memory location PLC 2
- F7: Delete
Delete the test program

7.15.2 Uploading Test Program from SD card (F2)


Fig. 7-34 | Program database: Load test program

With this function you can upload test programs previously archived on SD card back into the Sartocheck® 4 plus program database. Insert the SD card into the Sartocheck® 4 plus and press PAGE DOWN.


Fig. 7-35 | Program database: Load test program

The data on the SD card are uploaded and displayed in a program selection list. Here you can select the programs as described (see Fig. 7-28: Program databases: Internal memory).


Fig. 7-36 | Test program:
Assign memory location number

After you have selected the program or the program block, you will be prompted to assign a memory location number and | or memory location to this program or program block in Sartocheck® 4 plus.

Fig. 7-37 | Program database:
Program memory location

The test program and | or the test program section will be uploaded and saved to the selected internal program memory location.


Fig. 7-38 | Test program:
Assign memory location number

If the memory location number you have selected is already assigned, a warning will be displayed. With PAGE UP, you can select a new program memory location. By pressing Enter, the program data already stored in the internal memory will be overwritten with the new data.

Loading the test program from the network (F3) follows the same procedure as loading from an SD card (F2).

7.15.3 Saving Backup to SD Card (F4)


Fig. 7-39 | Program database: Save backup

With the backup function, all test programs in the internal memory will be written on an SD Card (F3) and/or the test programs in the internal memory will be rewritten (F4). Insert the SD Card into the SD card reader and then press PAGE DOWN.

Only use the backup function for saving data, not for archiving data. The backup files can only be viewed with Sartochek® 4 plus because file compression is specific to Sartochek® 4 plus. For archiving your program files, we recommend you use the (F1) function "Test programs from internal memory" -> (F3) "Save on an SD card."


Fig. 7-40 | Program database: Save backup

You have to individually name your backup files. Type in the backup file name and confirm by pressing enter.


Fig. 7-41 | Program database: Save backup

Then all data in the program database will be written to the backup file.

If the backup SD card already contains a backup file, you must confirm that the data will be overwritten before saving or cancel the backup routine.

Saving the backup to the network (F5) follows the same procedure as from an SD card (F4).

7.15.4 Loading Backup from SD card (F6)

When reloading backup files, insert the backup SD card into the SD card reader.


Fig. 7-42 | Program database: Read backup


Fig. 7-43 | Program database: Load backup

If you want to save more than one backup file, open the backup selection list, where you can select the file you want.


Fig. 7-44 | Program database: Load backup

When re-saving a backup file, all test programs in the internal memory are overwritten with the backup files.

Loading the backup from the network (F7) follows the same procedure as from an SD Card (F6).

7.15.5 Deleting PLC Programs (F8)


Fig. 7-45 | PLC program deletion

This function will delete programs in memory locations PLC 1 and PLC 2 from the selected test unit.

7.16 Visualize (F3 – Main Menu)


Fig. 7-46 | Visualize

This function will display the test status of any currently active test unit. Only active test units will be selectable. Direct access to this screen is also possible by using the status buttons at the bottom of all screens.

7.17 Results Database (F5 – Main Menu)


Fig. 7-47 | Results database

The results database is structured analogously to the program database. This is where all test results are saved. The functions of the results database are similar to those of the program database.

- F1: Test results from internal memory
This will display all test results inside the Sartocheck® 4 plus Filter tester
- F2: Write test results from SD Card to memory
This enables the results to be uploaded to the internal memory of the Sartocheck® 4 plus Filter tester
- F3: Upload test results from network
Displays all test results on the network in the corresponding directory
- F4: Save all test results to SD card
Saves all test results in the Sartocheck® 4 plus Filter tester on an SD card
- F5: Save all test results to network
Saves all test results on the network in the corresponding FTP directory
- F6: Manage audit trail
Manages the audit trail and defines which events are traced

With F1, all test results from the internal memory are displayed in a selection list. The list is structured as follows:

No.	Date	Serial no.	Type Name	Res
1	22.11.2005 12:12	17602314	DIF DIF	NN
2	22.11.2005 9:00	17602314	DIF DIF	OK
3	22.11.2005 9:02	17602314	DIF DIF	NN
4	14.11.2005 17:53	17602314	DBP Komplett	OK
5	14.11.2005 17:58	14801201	DBP Komplett	NN
6	14.11.2005 17:53	17602315	DBP Komplett	NN
7	14.11.2005 17:58	14801201	DBP Komplett	NN
8	14.11.2005 18:48	17602315	DBP Komplett	OK
9	14.11.2005 18:48	17602314	DBP Komplett	NN

Fig. 7-48 | Results database: Internal test results

The first column shows whether the result has already been saved. The saved program displays a diskette symbol. In the following columns, the current memory location number, the date of creation with time, serial number of your Sartocheck® 4 plus Filter tester, the integrity test type as well as the individual name of the test program are displayed. The last column provides information about the test evaluation:

- OK: Test passed
- NN: Test cancelled, no result
- F: Test not passed, line is displayed in red.

You can highlight the results by entering the list number or using the touch screen. Above the results list there is an entry line "From ### to ###". In this line, you can highlight the test results block.

Confirm your selection by pressing enter. The highlighting is cancelled by selecting "Remove highlighting."

No.	Date	Serial no.	Type Name	Res
3	22.11.2005 11:47	17602315	DIF DIF	NN

Fig. 7-49 | Results database: Internal test results 1

Above the function keys F1-F5, the highlighted test result is displayed. In the highlighted test results blocks, you can scroll through the selection list with LIST UP|DOWN. The following functions can be selected:

- F1: View
The test results are displayed
- F2: Printout
Prints out the test results on the internal printer
- F3: Save to SD card
Saves the test results to SD card
- F4: Save to network
Saves the test results to the network
- F5: Delete
Deletes the test results

In the menu Results database: Save to SD card and save to network

The highlighted result will be saved to the SD card or the network.

When saving the results, Sartocheck® 4 plus Filter tester establishes a file name. This file name is uniquely assigned. It is comprised of the following:

yymmdd_hhmm_04WIT12345678#.s4r

- yymmdd: Year, month, day in which the test was started
- hhmm: Hour, minute at which the test was started
- 04: Memory location number of the test program
- WIT: Test type
- 12345678: Serial number of the Sartocheck® 4 plus Filter tester
- #: Counting number (It is theoretically possible that a test is started twice in one minute, if the first test was cancelled immediately after it was started)

Fig. 7-50 | Results database: Save

When loading the test results from SD card or the network, all test results located on the SD card or network are displayed in a list. This list has the same structure as the list of the internal test results. Here you can highlight individual results or result blocks and display or print them out.

It is not possible to save external result files back into the Sartocheck® 4 plus Filter tester results database. It is only possible to display the external result files.

7.18 Special Functions


Fig. 7-51 | Special functions

In the menu "Special functions", you will find Sartocheck® 4 plus Filter tester functions that are (F6 – Main Menu) not directly linked to bag or filter testing. Here, the function tests for the management unit and the test unit are carried out. The cleaning procedure and drying procedure for the internal pneumatic block and internal volume reference tank can be selected here. Furthermore, pressure gauge functions can be activated in the special functions menu. The user can change their own password as long as they have access rights to this menu. For more details, refer to the menu "Service user management".

In the menu "Special functions", the corresponding functions can be activated with the F1-F6 keys.

While the Special Functions menu is active, visualizing a currently running test is not possible. Any test unit that is currently active will not be selectable through the Special Functions menu.

7.18.1 Function Test MU (F1)

With the Function test MU (management unit), you can check Sartocheck® 4 plus Filter tester peripherals.

Fig. 7-52 | Function test MU

- F1: Printer function test
Prints out all print characters in black and red on the internal printer
- F2: SD card function test (SD card reader)
After inserting an SD card, you can select from three different test procedures: Read only | Check file format | Write and read
- F3: Touch screen function test
The function of each individual touch field can be checked
- F5: Format SD card

7.18.2 Function Test TU (F2)


Fig. 7-53 | Special functions: Function test TU

The TU Function test will check the pneumatic block of the selected test unit. Here you are interactively guided through a test menu in which all valves and outlet ports are individually checked. Please follow instructions in the Sartocheck® 4 plus Filter tester display. Using the valve diagram you will be able to identify immediately which valves are working properly.

The indicator lights next to the valve will light up green. In the event of an error, they will light up red.

7.18.3 Cleaning (F3)


Fig. 7-54 | Cleaning

To be able to connect the cleaning kit to the Sartocheck® 4 plus Filter tester, please connect the test tubing between the Sartocheck and the cleaning kit appropriately. The cleaning menu allows you to flush all internal pneumatic parts completely, from the valves to the internal reference tank. The cleaning procedure is divided into multiple interactive steps. Sartocheck® 4 plus will guide you through each individual step. A diagram will show you which connections are to be connected with the flushing tubing during the current step. After you have installed all the connections according to instructions, press the start button. Sartocheck® 4 plus Filter tester then automatically carries out the current cleaning step. At the end of cleaning, the drying menu is activated automatically. We recommend that you also carry out a drying routine after every cleaning. This is done to ensure that no detergent solution remains in the pneumatic block.

If drying is not performed after cleaning, a message stating that the drying cycle has been cancelled will display every time a test is started. Once the drying cycle has been completed, this message will no longer display.


Fig. 7-55 | Cleaning

Please read the separate manual for the Sartocheck® Cleaning Kit carefully before starting the cleaning procedure! Only original Sartocheck® equipment may be used for cleaning and drying. Any misuse or use of unsuitable components may result in injury the operator or damage to the unit.


Fig. 7-56 | Cleaning



Fig. 7-57 | Cleaning

Description of the connections:

- a: IN – Inlet pressure tubing
Connect inlet pressure tubing
- b: Venting 2
Connect this outlet as defined
- c: OUT – Outlet pressure tubing to the filter housing
Connect this outlet as defined
- d: Venting 1
Connect this outlet as defined
- e: Ext.Ref.Tank / Venting 3
Connect this outlet as defined.
For the Sartochek 4 plus Filter tester you will need a coded Stäubli connector, see accessories.
- f: Shut-off valve
This valve may only be opened once all the cleaning equipment is installed and the cleaning procedure is started.
- g: Pressure Vessel
Connect the pressure vessel as defined and fill with cleaning fluid.
- h: Collecting Flask
Connect the collecting flask as defined. This vessel serves to collect the cleaning | rinsing fluid (drain) used.

The cleaning fluid we recommend:

Water (cold, warm or hot; never use warm or hot water for protein contamination!)
Buffer solutions, e.g. 0.9 % NaCl Solution
Sodium hydroxide solution, max. 1 M

In step 5, clean all the pneumatic components, i.e. the pneumatic block and the internal reference tank are flushed with water. This is done to ensure that no detergents remain in Sartochek® 4 plus Filter tester.


Fig. 7-57b | Cleaning

If you have not installed the connections in accordance with the diagram and the procedure has already started, we recommend that you cancel cleaning and restart. The pneumatic components can be damaged or their function impaired by improper connection. Improperly connected tubing can also impair the effectiveness of the cleaning procedure. Therefore, cleaning should be restarted.

Fig. 7-57c | Cleaning

The duration of the individual cleaning steps is predefined at 10 min per step. If it is necessary to change the time for the cleaning and drying steps, this can be done within the Service Menu.

To access this menu from the main menu, press F7 Service menu, F2 Software configuration, F10 Test unit parameters, F3 Edit parameters, F1 Cleaning | Drying.

7.18.4 Drying (F4)


Fig. 7-58 | Drying

After cleaning, drying follows. The drying program is structured analogously to the cleaning program. Here too, you will be interactively guided through the individual steps. The diagram will show you how the connections should be produced.

Description of the connections:

- a: IN – Inlet pressure tubing
Connect inlet pressure tubing
- b: Venting 2
Connect this outlet as defined
- c: OUT – Outlet pressure tubing to the filter housing
Connect this outlet as defined
- d: Venting 1
Connect this outlet as defined
- e: Ext.Ref.Tank | Venting 3
Connect this outlet as defined.
You will need a coded Stäubli connector, see “Accessories.”
- f: This valve must not be opened before the complete drying equipment is installed and the drying program started.
- g: Pressure Vessel
The pressure vessel is not used during the drying process.
- h: Collecting Flask
Connect the collecting flask as defined. This vessel is used for collecting the used cleaning | rinsing fluid.


Fig. 7-58a | Drying

7.18.5 Pressure Gauge Function (F5)

With the pressure gauge function, the internal and the external pressure sensor can be used for gas measurement. For measurement with the internal pressure sensor, the sensor (pressure gauge) is connected to the outlet. The remaining pneumatic components are turned off. Now, Sartocheck® 4 plus Filter tester works like a digital pressure gauge of accuracy class 0.1 and can be used to calibrate other pressure gauges.


Fig. 7-59 | Special function: Pressure gauge function

Do not use the pressure gauge function to monitor the pressure during liquid filtration!

During measurement with the internal pressure sensor, a connection needs to be made to the outlet connector to accurately read the atmospheric pressure.

7.18.6 Modify Own Password

Fig. 7-60 | Special Function:
Modify own password

According to the settings for the password expiration date, the user must change his or her password. By entering an old password, the user can store a new password and confirm. This password will then be valid until the defined expiration date is reached.

For more details, please refer to the menu:

Service -> Software configuration -> Management unit parameters -> Parameter processing-> Password

7.19 Service (F7 – Main Menu)


Fig. 7-61 | Service Menu

In the Service Menu, the hardware and software configurations are made, calibration of the internal pressure sensor and external pressure sensor are carried out and the user of Sartocheck® 4 plus Filter tester and/or new users are managed.

- F1: Hardware configuration
 - Interface: Define interface and external connections
- F2: Software configuration
 - Language
 - Date format
 - Time
 - Reference pressure data
 - Unit of measure
 - Record data
 - Management unit parameters
- F3: Calibration and adjustment
 - Internal pressure sensors and external pressure sensor inlets
- F4: User data management
 - Assign new user
 - Delete user
- F6: Manage user file
 - This function makes it possible to save and upload user rights to and from the net work and SD cards.

While the Service menu is active; visualizing a currently running test is not possible. Any test unit that is currently active will not be selectable through the Service menu.

7.19.1 Hardware Configuration
(F1-Service)

- F1: External reference volume
- F2: Network settings
- F3: Test unit address allocation
- F5: PLC interface test


Fig. 7-62 | Service Menu: Hardware configuration

7.19.1.1 External reference volume

In the hardware configuration, the volume for the external reference tank for each test unit is set separately.

The external reference tank function is only valid for the Filter tester.


Fig. 7-63 | Service Menu: Hardware configuration

7.19.1.2 Network settings


Fig. 7-64 | Network setting

Fig. 7-65 | Network connection

Selecting F1 from the Network settings menu brings up the Network connection menu. Within this menu is where the IP address of the Sartocheck® 4 plus Filter tester is identified as well as the addresses of the server and gateway. Explanations of the parameters follow:

F1: Address allocation

There are two (2) choices for this item

Automatic: This is chosen if a dynamic address assignment via DHCP is made. If this option is chosen, all subsequent fields will be grayed out, as the DHCP server will automatically provide the parameters.

Manual: If a static address assignment is being made, then this option should be chosen. At in minimum the SC4 plus IP address and the subnet mask must be entered.

F2: SC4 plus IP address

If a static address assignment is selected (manual), the IP address of the Sartocheck® 4 plus Filter tester is entered here

F3: Subnet mask

Use the information from the subnetwork screen to define the network type

F4: DNS server address

If a static address assignment is selected (manual), enter the IP address of the DNS server potentially existing on the network.

This is required for the resolution of addresses such as ftp.company.com

F5: Gateway

A gateway address is needed if the FTP server is located in a network other than the Sartocheck® 4 plus Filter tester.

Fig. 7-66 | FTP Server

Selecting F2 from the Network settings menu brings up the FTP server menu.

Within this menu is where the IP address and port of the FTP server are set, as well as the user name and password that the Sartocheck® 4 plus Filter tester should use when accessing the FTP server. Explanations of the parameters follow:

F1: Server address

Enter the FTP server address as the numerical IP address (xxx.xxx.xxx.xxx) or, if a DNS server is available, as a name such as ftp.customernet.com

F2: Port

Enter the FTP port (typically port 21) to which the request should be made for the connection setup

F3: User name

Enter the user name of the user setup on the FTP server

F4: Password

Enter the applicable password

F5: Timeout

Enter the time out (or wait time) for responses from the FTP server


Fig. 7-67 | FTP directories

Selecting F3 from the Network settings menu brings up the FTP directories menu. Entered here are the names of the directories that have been created in the user directory on the FTP server. Sub directories can also be entered, separated by the slash character (results/today)


Fig. 7-68 | Network settings

Selecting F4 from the Network settings menu will allow the user to test the configuration of the network connection.

Fig. 7-69 | Periodic save

Selecting F5 from the Network settings menu brings up the periodic save menu. This menu is used to activate and configure the periodic save feature. Explanations of the parameters follow:

- F1: Periodic save
Select "Active" or "Not active" to turn on or off this feature
- F2: Time Period
Select a predetermined time period or manual
- F3: Manual time input
If manual was selected above, enter a time between 5–10000 minutes

7.19.1.3 Test unit address


Fig. 7-70 | Service Menu: Hardware configuration

This function sets the address of the test unit so that the Sartochek® 4 plus Filter tester can allocation recognize it. The address Test Unit #1 is reserved for the Sartochek® 4 plus Filter tester.

Allowable addresses for the MultiUnits are Test Unit #2 through Test Unit #5. MultiUnits come from the factory with no address allocation, and therefore must be set by the user.

When connecting MultiUnits to the Sartochek® 4 plus Filter tester for the first time, it is important that only one new MultiUnit be connected at a time. Two or more units that have the same address allocation, or no address allocation (default setting from the factory), will cause a communication conflict.

When transferring a MultiUnit to be used with a different Sartochek® 4 plus Filter tester, it is advisable to first set the address allocation back to the default settings.


Fig. 7-71 | Service Menu: Hardware configuration

MultiUnit address allocation procedure:

- ▶ MultiUnit should be turned on, but the communication cable to the MultiUnit should not be connected.
- ▶ From the Hardware configuration screen of the Sartochek® 4 plus Filter tester, press F3 (Test Unit address allocation).
- ▶ Connect the MultiUnit communication cable when prompted by the Sartochek® 4 plus Filter tester and press "Page down".
- ▶ The new MultiUnit will be displayed in F7 as "New test unit" with the appropriate serial number.
- ▶ Select F7 and press "Page down".
- ▶ Available test unit address locations will be shown.
- ▶ Select the desired address location and press "Page down".

Change a MultiUnits address allocation:

- ▶ From the Hardware configuration screen of the Sartochek® 4 plus Filter tester, press F3 (Test Unit address allocation), then press "Page down".
- ▶ Select the appropriate test unit and press "Page down".
- ▶ Select the desired address location and press "Page down".

To set the MultiUnit back to the factory default setting, select F7 [Reset to standard address] and press "Page down".

7.19.2 Software Configuration
(F2-Service)


Fig. 7-72 | Service Menu: Software configuration

With the software configuration you can define:

Language	German English French Spanish Italian
Date format	dd.mm.yyyy dd/mm/yyyy dd-mm-yyyy mm/dd/yyyy yyyy/mm/dd yyyy.mm.dd
Date	as defined above
Time	hh.mm.ss
Reference	1000 mbar (abs)
Pressure	1013.25 mbar (abs)


Fig. 7-73 | Service Menu: Software configuration

Units of measure	
Pressure	mbar psi hPa
WIT unit	ml/10min ml/min
Temperature	°C °F K

All units of measure can be reset to Sartorius Stedim Biotech values using F1.

Pressure:	mbar
WIT unit:	ml/10min
Temperature:	°C

Fig. 7-74 | Service Menu: Software configuration

Record data
For each type of test (Diffusion, Bubble Point, etc.) you can define your basic record data. These defined parameters will then apply to all test programs. Go from the Service Menu – Software Configuration – Record data to the overview screen of the test type. Now select the desired test type, e.g. Diffusion or Bubble Point Test. You are immediately taken to the list of the record fields.

Starting with company, building and department, you can now define individual record fields by placing a checkmark or cancelling in the “Selection” field.

Fig. 7-74b | Service Menu: Software configuration
(continued)

Fig. 7-75 | Service Menu: Software configuration

Fig. 7-76 | Service Menu: Software configuration


Fig. 7-77 | Service Menu: Software configuration

The selections of the individual fields determine which fields are displayed and accessible during an integrity test, as well as which fields will be printed out on the Sartochek® 4 plus Filter tester's internal printer.

Note: Fig. 7-75 is valid for the protocol settings of the combined diffusion and bubble point test.

Management Unit Parameters

In the submenu "Management Unit Parameters", you can set the following options:

- F3: Edit parameters
Edits the management unit software parameters
- F4: Save parameters to SD card
Saves the management software configuration file on SD card
- F5: Save parameters to network
Saves the management software configuration file on the network in the corresponding FTP directory
- F6: Load parameters from SD card
Loads the management software configuration from SD card
- F7: Load parameters from network
Loads the management software configuration from the corresponding FTP directory in the network
- F8: Print parameters
Prints the parameters of the management unit software
- F9: Set standard parameters
Sets back the factory default parameters


Fig. 7-78 | Service Menu: Software configuration

Edit Parameters

F1: Password

You can Activate | Deactivate the password expiration date. When the password expiration date is activated, the password must be renewed after the password expiration date. Minimum duration is 1 day and a maximum 366 days. After the password expiration date has expired, you will automatically be prompted to renew the password. The max wrong password attempts can be set from 1-999. The minimum password length can be set from 1-999. Change password new user: if activated, the newly created user must change the password during the initial connection.

F2: Test program

Automatic Test Time Mode can shorten the test time if the values are constant within a defined tolerance. When Automatic Test Time Mode is deactivated, the test time is used that is set during the programming of the test program. The sequential tests function is used for repeating the same test.

F3: Hardware

Furthermore, the internal printer can be turned ON | OFF.

Management unit parameters

– Save parameters to SD card, F4

– Load parameters from SD card, F5

With these two options, you can save Sartochek® 4 plus management unit parameters to an SD card, and | or load saved data from an SD card into Sartochek® 4 plus Filter tester. The management unit parameters stored in Sartochek® 4 plus Filter tester are overwritten with those on the SD card.

F4: Auto Log Off

To activate | deactivate the auto log off and to define the timeout

F5: Results database

To define after how many test results the Sartochek® 4 plus Filter tester will give a warning message that the memory will soon be full

F6: PDF

To activate | deactivate saving of the test result files in PDF format.

To select the PDF page format A4 | Letter

F7: Data input

Input device: Keyboard/Keyboard+barcode/Only barcode scanner

Report field: select which field should be scanned

for recognition by the barcode


Fig. 7-79 | Service Menu: Software configuration

Here you can activate and deactivate the Auto Log Off function and set the auto log off time between 1 and 99 minutes

Fig. 7-80 | Service Menu: Software configuration

Defines when the Sartocheck® 4 plus Filter tester gives a warning message that the warning limit has been reached. Saving the test results on an SD card or on the network followed by deleting them is then recommended.


Fig. 7-81 | Service Menu: Software configuration

When uploading from the network or from an SD card (F6 or F7 from "Management unit parameters" menu), you get a warning message to prevent involuntary overwriting of data


Abb. 7-81a | Service Menu: Software configuration

Save as PDF file

When data are stored on an SD card or via a network, the option of saving PDF files can be activated and the format (A4 | letter) can be defined. The content of the PDF file is identical to that on the paper printout.

7.19.3 Calibration | Adjustment
(F3-Service)

Calibration serves to compare the internal pressure sensor and the external pressure sensor with an appropriate reference. If the deviations are outside the defined tolerance range, an adjustment is required. During adjustment to the defined pressure and | or current, surge sites are set and compared with the reference. The reference values are entered into the adjustment chart. After the adjustment is completed, you can immediately start a new calibration.

The calibration curves are stored in Sartocheck® 4 plus Filter tester and can be printed out at any time.

For the internal pressure sensor and the external pressure sensor, you set a calibration interval and a calibration reminder.

Before performing the actual calibration, allow Sartocheck® 4 plus Filter tester to warm up to operating temperature. This is done by turning Sartocheck® 4 plus Filter tester on 30 minutes prior to calibration.

The procedures for calibration and adjustment are identical. The calculation procedures differ in that calibration calculates and displays the deviation to internal adjustment curves. During adjustment, on the other hand, a new adjustment curve is plotted.

Firstly, you define which sensor and | or sensor inlet you want to calibrate. In the case of multi-unit operation, all logged-on multi-units are displayed. Their sensors are also calibrated via the central unit and adjusted, if applicable. By contrast, all calibration curves are stored in each multi-unit.

Use the function keys F1-F5 to select the test unit (MultiUnit). The internal pressure sensor calibration values of the respective test unit may be printed.

Use the "Print" button to print out the calibration protocol of the sensors highlighted in blue.

Press "Page Down" to go to the next window, in which the calibration interval and the calibration reminder for the sensor selected can be predefined.

The calibration should be checked at least once or twice annually. For your assistance, Sartocheck® 4 plus Bag tester provides the option to enter the calibration interval and a calibration reminder. These entries only serve as a reminder and do not affect the functioning of the device.

The calibration interval can be preset from 1 day to a maximum of 366 days. The same applies to the calibration reminder, which can be predefined from 0 to 366 days. The value "0" disables the reminder function.

After you have set the calibration interval and the reminder interval before the calibration, press Page Down to go to the next screen page. Use F1 to define whether a calibration or adjustment is to be performed. You can also enter information about the reference equipment (F4-F10 for internal pressure sensor, F4-F8 for external pressure sensor).


Fig. 7-82 | Calibration|Adjustment


Fig. 7-83 | Calibration

7.19.3.1 Calibration and Adjustment of the Internal Pressure Sensor

Sartocheck 4plus V. V2001 SN:12345678 User name: Mentor

Service menu: Calibration/adjustment

TU5: Internal pressure sensor

F1 Function: Calibration

F2 Range: 0-6 bar

Reference sensor

F5 Type: Absolute pressure sensor

F6 Model name: ABC

F7 Serial no: 123456

F8 Certificate no: 0087654

F9 Internal no: 54321

F10 Comments: none

Main Menu

Page up

Page down

Enter

Test unit 1: Modify

Test unit 2: Off

Test unit 3: Off

Test unit 4: Off

Test unit 5: Modify

Fig. 7-84 | Relative pressure reference

Connect the Sartocheck 4 plus Filter tester to a supply with sufficient current and compressed air (max. 7–10 bar). Connect, the outlet "OUT" to your reference sensor. Enter the data of your reference in the designated list. In the event that you have a relative pressure sensor as reference, enter the data in the list on the next page.

Service Menu Calibration | Adjustment

- TU1: Internal pressure sensor
- F1: Select operational step adjustment
- F2: Enter the following pressure range calibration:
- Bag tester only: 0 – 0.6 bar
 - Filter tester: 0 – 6 bar or 0 – 8 bar

Reference sensor

- F5: Type of sensor (not selectable for the Bag tester)
- F6: Identification (e.g. manufacturer)
- F7: Serial number
- F8: Certificate number
- F9: Internal number
- F10: Remarks

The information for describing the reference sensor has no effect on the actual calibration | adjustment, but is entered for documentation purposes and is printed out on the records accordingly.

Sartocheck 4plus V. V2001 SN:12345678 User name: Mentor

Service menu: Calibration

TU5: Internal pressure sensor

F1 Ambient pressure: 350 mbar(abs)

Ambient pressure reference

F4 Model name: Reference-P

F5 Serial no: 123456

F6 Certificate no: 0087654

F7 Internal no: 0123456

F8 Comments: none

Main Menu

Page up

Page down

Enter

Test unit 1: Modify

Test unit 2: Off

Test unit 3: Off

Test unit 4: Off

Test unit 5: Modify

Fig. 7-85 | Ambient pressure reference

If a relative pressure sensor is used as reference (only valid for the Filter tester), enter the additional specifications for the ambient pressure reference into the corresponding list.

Service Menu calibration

- TU1: Internal pressure sensor
- F1: Unit of measure

Ambient pressure in mbar (abs); only shown if the relative pressure sensor has been selected

- F4: Identification of reference current source (e.g. manufacturer)
- F5: Serial number of reference
- F6: Calibration certificate number
- F7: Internal number
- F8: Remarks

With Page Down, start the actual calibration and | or adjustment.

Sartocheck 4plus V. V2001 SN:12345678 User name: Mentor

Service menu: Calibration

TU5: Internal pressure sensor

Nominal press.	Reference press.	mbar(abs)
1502	1502	mbar(abs)
2502	2502	mbar(abs)
3502	3501	mbar(abs)
4502	4502	mbar(abs)
5502	5502	mbar(abs)
6502	6502	mbar(abs)

Main Menu

Page up

Page down

Enter

Test unit 1: Modify

Test unit 2: Off

Test unit 3: Off

Test unit 4: Off

Test unit 5: Modify

Fig. 7-86 | Calibration

Sartocheck® 4 plus Filter tester now runs up to the first pressure level. After pressure stabilization, you are prompted to enter the pressure applied to the reference sensor into the field "reference pressure". Afterwards, the next pressure level is run automatically, and stabilized. You are again prompted to enter the reference pressure value. Continue to repeat this procedure until the last pressure level has been processed. Then, calculation of the data is performed and, depending on the procedure selected (calibration or adjustment) the results table is displayed.

The figure 7-82 shows the pressure steps for the Sartocheck® 4 plus Filter tester.


Fig. 7-87 | Result of the calibration

During calibration, the results are displayed as a deviation from the adjustment curves.

Sartochek® 4 plus Bag tester:

If the deviation is below 0.025 % FS, an adjustment is not required.
At deviations of 0.025 – 0.5 % FS, an adjustment is recommended.
At deviations > 0.5 % FS, an adjustment is mandatory.

Sartochek® 4 plus Filter tester:

If the deviation is below 0.05 % FS, an adjustment is not required.
At deviations of 0.05 – 0.1 % FS, an adjustment is recommended.
At deviations > 0.1 % FS, an adjustment is mandatory.


Fig. 7-88 | Result of the adjustment

During adjustment, a new adjustment curve is calculated and activated using the entered values. The data of the adjustment can be printed out on the internal printer after the adjustment is completed. The new adjustment curve can then be checked with a calibration.

“Page up” takes you back to the selection screen for test units (MultiUnits) and sensors. Depending on “F10 perform calibration | adjustment?”, “Page down” is used to immediately start an additional calibration or adjustment.

7.19.3.2 Calibration and Adjustment of the External Pressure Sensor Input

The calibration of the external pressure sensor input is only valid for the Sartochek® 4 plus Filter tester and the Filter tester MultiUnit. Connect the Sartochek® 4 plus to the current and compressed air supply. Connect the electrical outlet "Ext. Sensors" to the external reference current source (4-20 mA). Enter the data on your reference into the designated list (see Figure 7.34).

Fig. 7-89 | External reference

Service-Menu Calibration | Adjustment
F1: select Calibration | Adjustment

Current Reference

F4: Model name
F5: Serial No.
F6: Certificate No.
F7: Internal No
F8: Comments

With Page Down, start the actual calibration and | or adjustment.

Nominal value	Reference value	
4.000	4.000	mA
8.000	8.000	mA
12.000	12.000	mA
16.000	16.000	mA
20.000	20.000	mA

Fig. 7-90 | External pressure sensor input

Enter the corresponding target value of the reference current supplied and wait for approx. 3 s. If there is a difference between the predefined value and the value read by the reference, then enter the read-in value in the table. If the values match exactly, go with PAGE DOWN to the next reference point. Continue to repeat this procedure until the last reference point has been run up and processed. Calculation of the data is now complete. Depending on the procedure selected (calibration or adjustment), the results table is displayed.

Nominal value	Actual value	Reference value	
4.000	4.000	4.000	mA
8.000	8.000	8.000	mA
12.000	12.001	12.000	mA
16.000	16.001	16.000	mA
20.000	20.000	20.000	mA

Fig. 7-91 | Result calibration

During calibration, the results are displayed as a deviation of the adjustment curves.

This classified result is output in addition to the table providing the values for the current supplied. Use "Print" to generate a printout for the calibration just carried out.

Nominal value	Actual value	Reference value	
4.000	4.000	4.000	mA
8.000	8.000	8.000	mA
12.000	12.001	12.000	mA
16.000	16.001	16.000	mA
20.000	20.000	20.000	mA

Fig. 7-92 | Result adjustment

During adjustment, a new adjustment curve is calculated and activated using the entered values. The data of the adjustment can be printed out on the internal printer after the adjustment is completed. The new adjustment curve can then be checked with a calibration.

"Page up" takes you back to the selection screen for test units (MultiUnits) and sensors. "Page down" is used to immediately start an additional calibration or adjustment, depending on "F10 perform calibration | adjustment?".

7.19.4 User Data Management
(F4 – Service)

Fig. 7-93 | Service Menu: Manage user data

The menu item “Manage users” can only be accessed by users with administrator rights. Here, new users can be created or old users deleted. The user management is based on three user hierarchies:

- 1. Standard User
- 2. Administrator
- 3. Service (Sartorius Stedim Biotech)

You can define the user rights for 1. and 2. Hierarchy level 3 is only accessible to trained service technicians.

When assigning standard user rights, the System administrator can assign individual user rights from a “rights matrix”.

F5: Add new user
With F5 function, select “Add New User”.
Now enter the user ID. Keep your user ID as short as possible. Then enter the full user name and the user’s function.

For the login procedure, you need your user ID.

These fields are “must” fields. Here it is mandatory that entries be made.
When finished, enter the password of the new user. After repeating the password, press the “OK” button.
Now the new user is logged in.

You are automatically taken to the display “Manage users”.
In the line “F1 User name”, the name of the new user is displayed.
With PAGE DOWN, you scroll to the next screen.

F6: Delete user
This function deletes the user you have selected in F1 (User name). It is not possible to delete your own user name.

Fig. 7-94 | Service Menu: Modify user rights

F7: Modify user rights
After a new user has been added, you can now set up their special access rights. Using the table, select the rights you want with the F1-F8 keys. The checkmarks next to the rights in the list have been assigned to the user. Checkmarks are placed or deleted with F1 to F10 or directly by touching the touchscreen.

In case a user account has been blocked due to too many attempts, the number of attempts can be reset to zero with or without the need to change the password.

F8: Reset password count
Resets the previous number of attempts to zero.

F9: Enter new password
This field allows specifying the new password for a given user.

F10: Repeat new password
Use this field to repeat the new password entered in F9.

7.20 Using the Barcode Scanner (optional)


Fig. 7-95 | Activation of the barcode scanner

Programming a Test and Automatic Test Selection after Scanning

In order to make the process of data entry more convenient and error-free, a barcode scanner can be connected to Sartochek® 4 plus Filter tester. Please note that it is only allowed to use official Sartochek accessories (barcode scanner type 26288---BS) to operate the unit.

The barcode scanner can be connected to the port "Serial Port MU" on the backside of the unit. The scanner will be immediately recognized and is ready to use ("plug & play").

The function "Keyboard + barcode scanner" is activated by default upon delivery. It does not oblige the operator to use a barcode scanner. Data entry can also be done using the keyboard only.

Before using the barcode scanner, this function must be activated in the Service Menu:

Principally, the Sartochek® 4 plus Filter tester barcode scanner is able to read all 1D-barcodes and 2D-data matrix codes. During the process of programming a test, the scanner can be used to enter the data into the respective comment fields. For this purpose, it is recommended to create a sheet with all the required coded information (e.g. company name, name of the building etc.). This will increase the speed of creating test programs.

Use of the barcode scanner allows automatic program selection after scanning. This feature can be activated | deactivated separately (compare Fig. 7-91). This means that Sartochek® 4 plus searches within the program database for information which matches the information coded in the barcode. After scanning, the unit displays only those test program(s) suitable for the scanned filter cartridge. When numerous test programs are defined in the program database, the process of finding the suitable test is improved and mistakes are avoided.

Sartochek® 4 plus Filter tester allows you to either search within a single comment field or to use all fields to finding the correct information. (compare Fig. 7-91).

The following procedure is recommended:

1. Defining a test program in the program database: it is recommended to use the barcodes of the filter cartridges during test programming and to enter the coded information into the dedicated comment field by scanning the barcode. Please note that the cursor has to be in the correct comment field before the cartridge is scanned. Following this procedure, more than one test program for each cartridge type can also be programmed (e.g. separate programs for BP and Diff Test).
2. Starting a test: To start a suitable test for a specific filter cartridge, please choose "Test Start" (F1) from the main menu. Then scan the barcode of the filter cartridge (exactly the barcode which was used for programming the test). After the process of scanning, Sartochek® 4 plus Filter tester will only display the suitable test program(s). If only a single test program can be found, Sartochek® 4 plus Filter tester will automatically select that one.
3. Continue according to the standard way of starting a test.

8. Maintenance | Servicing

Maintenance Intervals

- Every 6-12 months

To ensure that your Sartocheck® 4 plus Filter tester will give you many years of high-quality, consistent service, you can have it regularly calibrated. This is a service we offer as part of a maintenance agreement. If the Sartocheck® 4 plus Filter tester should show any unusual inconsistencies during an integrity test sequence, please notify your local Sartorius Stedim Biotech office or directly contact the corporate headquarters in Germany:

Sartorius Stedim Biotech GmbH
Service | FT
August-Spindler-Strasse 11
37079 Goettingen, Germany
Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius-stedim.com

Whenever you need to contact our technical customer service or require spare parts, please contact your Sartorius Stedim Biotech representative or Sartorius Stedim Biotech GmbH directly with the name and identification of the unit.

Before exchanging the fine-wire fuses, turn off the unit and unplug it from the wall outlet (mains supply).

8.1.1 Exchanging the Fine-wire Fuses

Fuses:

Two 250 V/16 A time-delay fuses compliant with DIN 41662 IEC 60127-2/III or IEC 127-3.

Remove the fuse socket and pull out the two fuses. Install the new fuses back in the connection block. Now plug the unit back into the wall outlet.

8.1.2 Sartorius Stedim Biotech Service Menu

The Service Menu contains all internal operating functions, limits and test abort parameters. Major errors in measurement can occur if they are changed incorrectly. Therefore, only a Sartorius Stedim Biotech service technician or a specially trained member of your staff should change the parameters in the Service Menu.

Service Menu instructions are only given to specifically named and authorized Service technicians who must confirm their receipt.

8.2 Cleaning

For cleaning, use a wet cloth and wipe afterwards with a dry cloth. Make sure that no humidity or moisture enters the device. Do not use solvents!

Follow the shutdown procedure (see Chapter 4.4) before carrying out maintenance, cleaning and repair work!

9. Troubleshooting

The work and instructions for “Troubleshooting” described in these operating instructions are explained in such a way that they can be understood and carried out by trained workers and qualified technicians (see Definition in Chapter 3.3 “Safety Precautions”) in

- Electrical | Electronics
- Mechanical | Maintenance

9.1 Disposal and repair instructions

This personnel must be provided with the appropriate tools and testing devices. In particular-, when maintenance and repair work is being carried out, the shutdown procedures (see Chapter 4.4) must be carried out under all circumstances. If the measures listed do not succeed, please contact the Sartorius Stedim Biotech GmbH.

9.2 Troubleshooting

Whenever errors occur, they appear on the display. Proceed as prompted by the menu guidance and observe the following instructions:

The current procedure continues after confirmation. The user is prompted with the following instructions:

- The current procedure has been stopped. A user decision is required and must be confirmed before Sartocheck® 4 plus continues to carry out or cancels the procedure.
- The current procedure has been cancelled. The user receives the error code “Enter elimination note” or the prompt to inform Sartorius Stedim Biotech Service.

If it becomes necessary to contact Sartorius Stedim Biotech Service, please have the following information ready:

- Serial number
(see the back panel of Sartocheck® 4 plus Filter tester test print out with header display)
- Version
(see header display)
- The four-digit error code

Error	Cause	Solution
The Sartocheck® 4 plus Filter tester Filter tester cannot be turned on.	Main switch is in the "Off" position	Turn on main switch
	Safety fuses blown	Have a Sartorius Stedim Biotech Service technician change the fuses.
The test cannot be started.	Tests can be only started if the menu items "Start test", "Program test" and "Test repeat" are activated	
Error Message: Temperature of the unit too high Temperature of the unit too low	The internal temperature of the unit is outside the limits set in the Service Menu T _{min} : +15°C T _{max} : +35°C	Wait for the temperature to adjust
No printout is generated	The paper is jammed in the printer	Remove paper jam
	There is no printer paper	Insert new roll of paper
	Automatic printing is turned off	Turn on internal printer.
The measured data obtained for the diffusion rate, pressure drop or intrusion are unexpectedly high, and or drop curve is non-linear	The filter element has not been wetted or has been incompletely wetted	Completely rewet the filter element
	There is housing or system leakage	Check system for leakage and eliminate filter system leakage
	Excessive temperature drift	Avoid temperature drift
	The test pressure is too high (already in the B.P. range)	Adjust the correct test pressure
	Stabilization and or test phase too short	Check test times and extend them, if applicable
	The water used for the intrusion test does not have sufficient surface tension	Make sure that the surface tension of the water is not < 70 dyn/cm (for the water intrusion test).
The bubble point is unexpectedly low.	The filter element has not been or has been incompletely wetted	Completely rewet the filter element
	There is housing or system leakage	Check system for leakage and eliminate filter system leak
	The wrong test code was entered for an individual integrity test	Enter the correct test code
	The surface tension of the wetting fluid was too low	Reset test parameters if you use a non-standard wetting fluid or contact the Validation Service of the Sartorius Stedim Biotech GmbH. Check temperature of wetting fluid.

Error-List

Legend for 'Type'

'R' Commit severity 'Remark'

'W' Commit severity 'Warning'

'E' Commit severity 'Error'

Placeholder in 'Usertext'

<ARG1>, <ARG2>, <ARG3> are placeholder that are filled with variable information like filename, testprogramnumber, testunit-number etc.

Err no.:	Type:	Usertext:
-121	E	File '<ARG1>' has incorrect checksum!
-120	E	The user management file is corrupt. Standard parameters will be set. Please change (and add) parameters in the service menu and save them!
-119	E	The contents of the file does not match the expected file structure.
-118	E	The stored calibration date of the external pressure sensor inlet is invalid
-117	E	The stored calibration date of the internal pressure sensor is invalid
-116	E	Error when writing the calibration date of the external pressure sensor inlet
-115	E	Error when reading the calibration date of the external pressure sensor inlet
-114	E	Error when writing the calibration date of the internal pressure sensor
-113	E	Error when reading the calibration date of the internal pressure sensor
-112	E	Audit trail file is damaged. Incorrect checksum. Please archive data to create a new audit trail file.
-111	E	Time exceeded when stopping the test unit
-110	E	The system log file could not be stored!
-109	E	The configuration file of the test unit is corrupt. The standard parameters will be set. Please change parameters in the service menu and save them!
-108	E	The configuration file of the management unit is corrupt. Standard parameters will be set. Please change parameters in the service menu and save them!
-107	E	The stored calibration protocol of the external pressure sensor is damaged. Please calibrate the external pressure sensor.
-106	E	The stored calibration protocol of the internal pressure sensor is damaged. Please calibrate the internal pressure sensor.
-105	E	The test program '<ARG1>' could not be printed. Error: '<ARG2>'
-104	E	The test program '<ARG1>' could not be deleted.
-103	E	The test program '<ARG1>' could not be saved. The SD card is full.
-102	E	The test program '<ARG1>' could not be saved.
-101	E	The test result '<ARG1>' could not be printed. Error: '<ARG2>'
-100	E	The test result could not be deleted.
-99	E	The test result could not be saved to the SD card. The SD card is full.
-98	E	The test result could not be saved.
-97	E	The log file could not be saved to the SD card. The SD card is full.
-95	E	Calibration interval of the external pressure sensor inlet invalid
-94	E	Calibration interval of the internal pressure sensor invalid
-93	E	Error when writing the calibration interval of the external pressure sensor inlet
-92	E	Error when reading the calibration interval of the external pressure sensor inlet
-91	E	Error when writing the calibration interval of the internal pressure sensor
-90	E	Error when reading the calibration interval of the internal pressure sensor
-89	E	Test results are damaged. The checksum is incorrect!
-88	E	Cannot close touch screen file! Error: <ARG1>
-87	E	Cannot open touch screen file! Error: <ARG1>
-86	E	Error when setting touch screen coordinates. Error: <ARG1>
-84	E	Error when stopping the test unit
-83 ... -82	E	Calibration data of the internal pressure sensor invalid

Err no.:	Type:	Usertext:
-81	E	Error when writing the calibration data from the external pressure sensor inlet
-80	E	Error when reading the calibration data of the external pressure sensor inlet
-79	E	Error while writing the calibration data from the internal pressure sensor
-78	E	Error when reading the calibration data from the internal pressure sensor
-76	E	New main memory cannot be allocated
-75	E	Unexpected file end when reading file
-74	E	Software versions of test unit, data management unit and communication server are incompatible Software version of test unit = <ARG1> Software version of data management unit = <ARG2> Software version of xbp2 server = <ARG3>
-73	E	System error file is too large
-71	E	Error when opening system error file
-70	E	Error when loading test result
-67	E	Program is damaged. The checksum is incorrect!
-66	E	Cannot terminate communication with the test units
-65	E	Cannot start communication with the test units
-64	E	System partition could only be mounted in readable condition!
-63	E	System partition could only be mounted in an unwritable condition!
-62	E	Cannot write to "sc4user" file. Please notify your system administrator
-61	E	Cannot read "sc4user" file. Please add users or notify your system administrator
-49	E	Cannot open audit trail file
-48	E	Cannot open error log file of the data management unit
-47	E	Cannot open error log file of the test unit
-46	E	Cannot open message window
-45	E	Unknown test method number '<ARG1>'
-44	E	Test unit number '<ARG1>' outside defined ranged
-43	E	Cannot open configuration file of test unit '<ARG1>' to write.
-42	E	Cannot open configuration file of test unit '<ARG1>' to read.
-41	E	Cannot open configuration file of processing unit to write to memory.
-40	E	Cannot open configuration file of processing unit to read from memory.
0 ... 2500	E	Internal error
6001	E	You cannot save any more test results. Please archive the result database.
6002	W	A total of '<ARG1>' out a max. '<ARG2>' results have been saved. Please archive the results database.
6003	E	There may be more than '<ARG1>' measuring points. A maximal of '<ARG2>' measuring points are allowed.
6004	E	Min. BP must be less than max. BP
6005	E	Limits for min. BP, max. BP: 100 mbar <= min. BP, max. BP <= 8000 mbar (1.45 psi < min. BP, max. BP <= 116.03 psi) mbar and hPa values have no decimal places, psi values have up to 2 decimal places.
6006	E	Limits for volume: 10 ml <= vol <= 150000 ml
6007	E	The test pressure must be less than the min. BP
6008	E	Limits for times: 1 min <= stability time, test time <= 99 min
6009	E	Limits for test pressure: 100 mbar <= test pressure <= 8000 mbar (1.45 psi < test pressure <= 116.03psi) mbar and hPa values have no decimal places, psi values have up to 2 decimal places.
6010	E	Limits for diffusion: 0.1 ml/min <= max. diffusion <= 999.9 ml/min Diffusion values can have up to 1 decimal place.
6011	E	Limits for pressure drop: 1 mbar <= max. pressure drop < test pressure (0.01psi < max. pressure drop < test pressure) mbar and hPa values have no decimal places, psi values have up to 2 decimal places.
6012 ... 6013	E	Limits for test water intrusion / water flow: 0.01 ml/min <= Max. water intrusion / max. water flow <= 99.99 ml/min (0.1 ml/10min <= Max. water intrusion / max. water flow <= 999.9 ml/10min) ml/min values have up to 2 decimal places, ml/10min values have up to 1 decimal place.
6014	E	Limits for min. and max. test pressure: 100 mbar <= min., max. test pressure <= 8000 mbar (1.45 psi < min., max. test pressure <= 116.03 psi) mbar and hPa values have no decimal places, psi values have up to 2 decimal places.
6015	E	Limits for MPD pressure increment: 25 mbar <= MPD pressure increment <= 500 mbar (0.36 psi < MPD pressure increment <= 7.25 psi) MPD pressure increment < 0.5*test pressure difference mbar and hPa values have no decimal places, psi values have up to 2 decimal places.

Err no.:	Type:	Usertext:
6016	E	Min. test pressure must be less than max. test pressure
6017	E	Limits for BP pressure increment: 25 mbar <= BP pressure increment <= 100 mbar (0.36 psi < BP pressure increment <= 1.45 psi) mbar and hPa values have no decimal places, psi values have up to 2 decimal places.
6018	E	Limits for BP FACTOR: 0.1 <= BP FACTOR <= 9.9
6019	E	Limits for brief times: 1 sec <= brief stabiliz. time, brief test time <= 999 sec
6020	E	Old password wrong! Please try to enter it again or notify your system administrator
6021	E	New passwords are not identical! Please try to enter them again or notify your system administrator
6022	E	New user ID already exists! Please try again
6023	E	New user name already exists! Please try again
6024	E	New password is same as old password! Please try again or inform your system administrator
6025	E	Your password has expired! Please enter a new password or inform your system administrator
6026	E	Wrong password! Please try again or inform your system administrator
6027	E	Unknown user! Please add a new user or inform your system administrator
6028	E	Log off not possible at the moment! Please try again later.
6029	E	Calibration of the internal pressure sensor is due. Please calibrate the sensor.
6030	E	Calibration of the internal pressure sensor is due in <ARG1> day. Please calibrate the sensor.
6031	E	Calibration of the external pressure sensor input is due. Please calibrate the sensor.
6032	E	Calibration of the external pressure sensor inlet is due in <ARG1> day. Please calibrate the sensor.
6033	E	The audit trail file is full! Please archive it!
6034	E	Limits for BP volume: 10 ml <= vol <= 50000 ml
6035	E	Error while setting the date/time for the data management unit
6036	E	Wrong parameter for setting the date
6037	E	Wrong parameter for setting the time
6038	E	Pressure differential after point '<ARG1>' is '<ARG2>' It must be at least as large as '<ARG3>'
6039	E	Diffusion differential after point '<ARG1>' is '<ARG2>' It must be at least as large as '<ARG3>'
6040	E	<ARG1> - <ARG2> error when creating program no.!
6041	E	Timeout while setting the linearization standard parameters
6042	E	Data not automatically stored '<ARG1>' : <ARG2>
6043	E	Serious internal error <ARG1> ; <ARG2> ; <ARG3>
6044	E	The test program '<ARG1>' could not be saved. <ARG2>
6045	E	Error accessing temporary directory
6046	E	The test result '<ARG1>' could not be saved. The SD card is full.
6047	E	Test unit address must be between range of 1 to 5.
6048	E	You do not have sufficient access rights to perform this operation!
6049	E	No user is logged in. Please log in first!
6050	E	No test has been run yet!
6051	E	Filter/system was vented.
6052	R	Test passed.
6053	R	Operation successful.
6054	R	Operation not successful.
6055	R	Test failed.
6056	E	No test result.
6057	E	Print failed, error <ARG1>.
6058	E	No test programs were found.
6059	E	The requested backup could not be loaded!
6060	E	The directory could not be transferred - Please check network connection.
6061	E	The temporary backup file could not be created.
6062	E	The backup file could not be transferred.
6063	E	The value entered is outside the valid range or number of decimals is wrong: <ARG1>

Err no.:	Type:	Usertext:
6064	E	The cleaning procedure was terminated early. - Close shutoff valve (f) - Disconnect "Venting 1" tubing from pressure tank (g) Please run drying procedure now.
6065	E	The drying process was stopped. Before the next test series, this process must be completely repeated. - Close shutoff valve (f) - Disconnect "Venting 1" tubing from pressure tank (g) - Disconnect remaining tubing from pressure tank (g). Caution: Tank (g) is being vented - Please disconnect all tubing from the connectors b-e.
6066	E	File number must be between 0 and 5!
6067	E	File '<ARG1>' cannot be transferred! '<ARG2>'
6068	E	There is no SD card inserted in the drive.
6069	E	File "<ARG1>" could not be copied!
6070	E	File "<ARG1>" does not exist!
6071	E	SD card test: User break
6072	E	SD card test: error!
6073	E	SD card test: error or no SD card in the drive.
6074	E	Formatting SD card: SD card is write-protected.
6075	E	Error occurred when formatting SD card
6076	E	Cannot establish network connection.
6077	E	Test result could not be transferred.
6078	E	Test result cannot be opened.
6079	R	The network connection functional.
6080	E	Cannot establish connection to server: <ARG1>
6081	E	Writing test unit timer equalization failed '<ARG1>'
6084	R	Test unit address allocation successful.
6085	E	Test unit address allocation not successful.
6086	E	The value entered is outside the valid range or number of decimals is wrong: <ARG1>
6087	E	Unexpected change in serial number. Expected: '<ARG1>' Now: '<ARG2>'
6088 ... 6089	E	Unexpected change in test unit status.
6090	E	Time exceeded when starting the program.
6091 ... 6092	W	File'<ARG1>' could not be transferred! '<ARG2>' Do you want to cancel the procedure?
6093	W	File'<ARG1>' could not be deleted! Do you want to cancel the procedure?
6094	W	Do you really want to change the measurement method?
6096	W	Do you really want to delete the entry selected?
6098	W	The current parameters are replaced by the factory settings. The parameters changed so far will be overwritten. Do you really want to do that?
6099	W	All data on the SD card will be overwritten. Are you sure you want to overwrite this data?
6100	W	This SD card is not suitable for saving Sartocheck4 files. Do you want to reformat the SD card?
6101	W	Drying is not performed completely! Do you want to start the test anyway?
6102	W	This cancels the ongoing procedure! Do you really want to do this?
6105	W	The file '<ARG1>' already exists. Do you want to overwrite it?
6106	W	There is no SD card in the drive! Please insert a formatted SD card in the drive to save the data. Do you want to repeat this procedure?
6107	W	There is no SD card in the drive. Please insert a formatted SD card to read the data. Do you want to repeat the procedure?
6108	W	There is no properly formatted SD card in the drive! Please insert a formatted SD card to save the data. Do you want to repeat the procedure?
6109	R	Test results not found.
6110	E	No function is currently assigned to the stop button. Either no tests are running or visualization is not activated.
6111	E	The maximum number of '%' users may not be exceeded.
6112	E	SD card is write-protected.
6114	E	It is not possible to start the selected procedure at the moment, because the former procedure has not yet been completed - Please try again later

Err no.:	Type:	Usertext:
6115	E	Not possible to start visualization, as no test is running
6116	E	Automatic program selection not possible as no program or report field input available
6117	W	-Cancel 'Error' and 'Ready' states; -Cancel 'Drying, internal/external adjustment not performed'; -Cancel 'Test repeat, PLC programs, last result'; -Abort program in progress Do you want this to be carried out?
6118	W	Once you start the following procedure, it must be performed completely. Do you want to start this procedure?
6119	E	Password was entered incorrectly too often. Please inform your system administrator.
6120	E	New user needs to enter new password.
6121	E	Temperature informations: <ARG1> <ARG2>
6122	E	<ARG1> ----- Temperature informations: <ARG2> <ARG3>
6123	E	No test result.
6124	E	Cleaning and drying record not available.
6125	E	Your password does not have the required minimum length of <ARG1> character(s); please enter new password.
6126	R	The drying process was completed. - Close shutoff valve (f) - Disconnect "Venting 1" tubing from pressure tank (g) - Disconnect remaining tubing from pressure tank (g). Caution: Tank (g) is being vented - Please disconnect all tubing from the connectors b-e.
6127	E	The user file of the management unit is corrupt.
6901 ... 6902	E	Printer server program cannot be addressed.
6903	E	Unknown print job.
6904	E	Job cannot be sent to printer server.
6905	E	Response cannot be recieved from printer server.
6906	E	Invalid print data.
6907	E	Print job list is full.
6908	E	Invalid test number for printing
6909	E	Cannot write print file.
6910	E	Cannot read print file.
6911	E	Insufficient main memory for sending print job
6912	E	Unknown error occurred during printing.
6999	E	Internal error! <ARG1>
8000 ... 8300	E	Communication error with test unit
9001	E	It is not possible to start the selected procedure at the moment, because the former procedure has not yet been completed - Please try again later
9004	E	Only the internal reference tank is active
9005	E	Both the internal and external reference tanks are active
9006	E	The basic configuration data were re-initialized: (Address of the test unit, serial number, software version of the test unit) - when the unit was initially operated at the factory - when the format of the data mentioned above was changed
9010 ... 9011	E	The raw data of the graphs were not sent to the processing unit because the required component currently does not exist
9012	E	The operating pressure ist too low. Please increase the operating pressure.
9013	E	The operating pressure ist too high. Please reduce the operating pressure.
9014	E	The operating pressure is not sufficient. Please increase the operating pressure.
9015	E	The operating pressure is not stable. Please wait.
9101	E	Error while venting Warning: The system is probably still pressurized!
9102	E	The required test pressure exceeds the calibration range of the pressure sensor - Recalibrate using a wider calibration range (0-8 bar/0-116.03 psi) - Work with a lower test pressure
9103	E	The +12 V power supply is out of tolerance.
9104	E	The -12V power supply is out of tolerance
9105	E	The +24 V power supply is out of tolerance
9106	E	Valve V1 does not close correctly.
9107	E	No compressed air supply is available or valve V1 does not open correctly

Err no.:	Type:	Usertext:
9108	E	The ambient pressure measured is out of tolerance
9109	E	The V2 valve for the reference tank is not operating properly or the reference tank is full
9110	E	The internal reference tank is full
9112	E	The test was subsequently cancelled because the operating pressure after the measurement was below 4 bar (58.02 psi)
9113	E	The volume of the valve unit is 0 ml
9114	E	During the pressure drop test, the pressure rose by more than 2 mbar (0.03 psi) - The temperature increased.
9115	E	During the pressure drop test, the pressure drop exceeded the limit - The system leaks
9116	E	The pressure drop test failed because the pressure drop exceeded the limit
9117	E	Internal error
9118	E	During calibration of the internal pressure sensor, the operating pressure for the selected pressure range is insufficient - Increase the operating pressure
9119	E	During calibration of the internal pressure sensor, all pressure values measured were the same so it is not possible to calibrate
9120	E	During calibration/adjustment of the internal pressure sensor, the difference between two reference pressures is too small for calibration to be performed.
9121	E	During pressure buildup, the nominal pressure was exceeded
9122	E	Timeout for pressure buildup exceeds the time limit
9123	E	Internal error
9124	E	During the stabilization time the system volume could not be stabilized sufficiently
9125	E	The pressure drop at the end of stabilization exceeds the limit
9126	E	Error while venting Warning: The system is probably still pressurized!
9127	E	The loss of pressure during volume measurement exceeds the limit - Please check system for leakage - Please check configuration
9128	E	The net volume determined is below 5 ml - Start cleaning cycle
9129	E	During volume measurement (with an external reference tank > 6000 ml), the net volume measured was so large that it was not possible to execute the test
9130 ... 9131	E	During volume measurement (without an external reference tank) the net volume measured was so large that it was not possible to execute the test
9132	E	The volume measurement yielded a net volume greater than 12 l for the bubble point test or 50 l for the complete test
9133	E	A considerable loss of pressure occurred at the beginning of the diffusion test - No filter housing was connected - The filter cartridge is missing
9134	E	During the diffusion test, a pressure increase occurred that may be caused by an increase in the temperature -Please wait for the temperature to equalize -Check the temperature of the wetting liquid
9135	E	During the diffusion test, the diffusion exceeded the limit -Please check the system for leakage -Check the configuration -Rewet the filter
9136	E	The last pressure value during the diffusion test exceeded the nominal pressure -Temperature increased during measurement
9137 ... 9138	E	During the bubble point or multi-point diffusion test, 5 consecutive pressure values exceeded the nominal pressure by 500 mbar (7.25 psi) -Reduce the operating pressure
9139	E	The diffusion at the first bubble point pressure increment exceeds the limit - Check the system for leakage - Enter the correct test class - Adjust the diffusion limit
9140	E	The diffusion at the second bubble point pressure increment exceeds the limit - Check the system for leakage - Enter the correct test class - Adjust the diffusion limit
9141	E	The diffusion from the third or higher bubble point pressure increment exceeds the limit - Therefore, the bubble point was not recognized - Check the system for leakage - Enter the correct test class - Adjust the diffusion limit
9142	E	Pressurization of housing has failed - Increase the operating pressure
9143	E	The internal power supply is partially defective
9144	E	The time it takes for stabilization exceeds the limit
9145	E	During calibration of the temperature sensor, the commands were sent in an invalid sequence
9146	E	During calibration of the temperature sensor, the type of sensor has been changed in between
9147	E	During calibration of the temperature sensor, two identical temperatures were measured
9148	E	During the bubble point test, the effective measuring time for the pressure curve was 0

Err no.:	Type:	Usertext:
9149	E	When testing for leakage during the water intrusion test, the water intrusion was too high: - No filter connected - No filter cartridge in filter housing
9150	E	The test program contains an invalid combination regarding the use of an external pressure sensor and an external valve - Please enter the correct combination
9151	E	During the water intrusion test, a pressure increase occurred that may be caused by an increase in temperature -Please wait for temperature compensation -Check the temperature of the wetting liquid
9152 ... 9154	E	Internal error
9155	E	The required stabilization time is too short
9156	E	During calibration/adjustment of the external pressure sensor input, the difference between two reference currents is too low so that calibration is not possible.
9157	E	During calibration of the external pressure sensor input, all pressure values/flow rates measured were identical so that calibration is not possible
9158	E	The time it takes to build up the pressure (to 100% nominal pressure) during external filling exceeds the limit
9159	E	At the end of external filling, the actual pressure in the filter housing deviates too much from the nominal pressure
9160	E	The duration of pressurization (up to 75% target pressure) during external filling exceeds the limit - The filter cartridge may be damaged - Connect pressure supply to external valves
9161	E	The duration of pressurization (up to 75% target pressure / filter fully wet)) during external filling exceeds the limit - Connect pressure supply to external valves
9162	E	During the bubble point and/or multipoint diffusion test, more than 200 pressure levels should be checked (e.g. [max. BP-min. BP] / pressure level > 200)
9163	E	Timeout for pressure drop exceeds the limit
9164	E	The external pressure sensor indicates too high a pressure value (possible short-circuit)
9165	E	Wrong code entered so access to service mode denied
9166	E	The required test was not activated; the unit is in the service mode
9167	E	The required service function was not activated; the unit is in the test mode
9168	E	The nominal pressure decreased during filling (75% nominal pressure > nominal pressure + water column) - Do not change pressure conversion vessel during filling - Reduce the "WIT filling factor" in the service menu
9169	E	The external pressure sensor has line defect or it is not connected
9170	E	The test has been discarded retrospectively as the external pressure sensor is not connected any longer
9171	E	During stabilization, nominal pressure exceeded target pressure
9172	E	During external filling, nominal pressure exceeded target pressure
9173	E	The flow measured is out of tolerance - select the required current and impress - if necessary, load factory-set adjustment first
9174	E	The ambient pressure measured is out of tolerance Warning: system pressure is probably still on!
9175	E	The test was cancelled by user/operator.
9176	E	The test was cancelled by the PLC.
9177	E	The test was cancelled by switching off the power supply.
9178	E	During the MPD test, diffusion exceeded limit A2 -Please check the system for leakage -Check configuration -Rewet filter
9179	E	Wrong parameters for the MDC test.
9180	E	The differential between external and internal pressure sensor is too large.
9181	E	The diffusion at the first bubble point pressure increment exceeds the limit - Check the system for leakage - Enter the correct test class - Adjust the diffusion limit
9182	E	The pressure drop in the bubble point pressure stage exceeds the limit - Check system for leakage - Enter the correct test class - Adjust the diffusion limit
9183 ... 9187	E	Internal error
9188	E	During volume measurement with an external reference tank, the net volume measured was so large that it was not possible to execute the test
9189	E	During volume measurement without an external reference tank, the net volume measured was so large that it was not possible to test the filter. It is recommended that you use an external reference tank.
9190	E	Volume measurement yielded a net volume greater than 150L
9191	W	Value below a defined min. limit - No filter/filter housing connected - Filter/filter housing incorrectly connected - Check configuration

Err no.:	Type:	Usertext:
9193	E	Value above the max. volume.
9194	E	Value below the min. volume.
9195	E	The function was cancelled.
9196	W	The internal temperature limit was exceeded.
9197	W	The environmental temperature fluctuations during Stabilization time/Test time are out of the permissible tolerance of $\pm$ <ARG1>. The actual temperature fluctuation is: Stabilization time: Min/Max (<ARG2>) Test time: Min/Max (<ARG3>)
9198	W	The internal and external temperature limits were exceeded.
9199	W	During monitoring using the pressure gauge function, the maximum allowed pressure was exceeded.
9201	E	An internal valve has short-circuited at the supply voltage
9202	E	An internal valve has short-circuited to ground
9203	E	An internal valve has a broken conductor
9204	E	An internal valve has short-circuited
9205	E	An external valve has short-circuited at the supply voltage
9206	E	An external valve has a broken conductor - Connect valve(s) - Electrically test valve(s)
9207	E	The external filling valve has a broken conductor - Connect valve(s) - Electrically test valve(s)
9208	E	The external drain valve has a broken conductor - Connect valve(s) - Electrically test valve(s)
9209	E	The external vent valve has a broken conductor - Connect valve(s) - Electrically test valve(s)
9210	E	The external block valve has a broken conductor - Connect valve(s) - Electrically test valve(s)
9211	E	An external valve has short-circuited
9212	E	An external proportional valve has short-circuited at the supply voltage
9213	E	An external proportional valve has a broken conductor - Connect valve(s) - Electrically test valve(s)
9214	E	The "small" external proportional valve has a broken conductor - Connect valve(s) - Electrically test valve(s)
9215	E	The "large" external proportional valve has a broken conductor - Connect valve(s) - Electrically test valve(s)
9216	E	An external proportional valve has short-circuited
9217	E	An internal proportional valve has short-circuited at the supply voltage
9218	E	An internal proportional valve has a broken conductor
9219	E	The "small" internal proportional valve has a broken conductor
9220	E	The "large" internal proportional valve has a broken conductor
9221	E	An internal proportional valve has short-circuited
9222	E	The external temperature sensor is not connected.
9301	E	Internal error
9302	E	internal error <ARG1>
9303	E	Internal error
9304	E	Valve has not been tested yet.
9305	E	Valve does not close.
9306	E	Valve does not open.
9307	E	Valve 3 does not close. Function test should not be continued.
9308	E	Valve 4 does not open. Function test should not be continued.
9309 ... 9310	E	internal error <ARG1>
9315	E	Communication error with A/D converter of the safety circuit
9316	E	Communication error for reference adjustment of the safety circuit
9320	E	The difference in pressure is too large. No meaningful adjustment possible.
9321	E	The safety valve does not open although it has been electrically deactivated.
9322	E	Adjustment could not be successfully completed within the time provided.
9323 ... 9325	E	The safety valve opened because the pressure has been exceeded.
9326	E	During test time, nominal pressure exceeded target pressure
9327	E	The safety valve circuit did not respond although the pressure was too high.

Err no.:	Type:	Usertext:
9328	E	The safety valve circuit was not adjusted.
9329	E	The safety valve circuit is not available.
9330 ... 9333	E	Communication error with A/D converter of the safety circuit
9401 ... 9406	E	Software versions of data management and test unit are incompatible
9407	E	Internal CRC error in test program
9408	E	Internal CRC error in test log data
9409	E	Internal CRC error in test result data
9410	E	Test results data not available.
9411	E	Software versions of data management and test unit are incompatible
9901 ... 9912	E	Internal error
9913	E	Date/Time invalid
9914 ... 9980	E	Internal error
9981	E	The test method selected cannot be performed on the test unit chosen.
9982	E	The selected function cannot be executed on the selected test unit.
9983 ... 9982	E	Internal error

10. Emergencies

The unit or any of its parts must be removed and disposed of in accordance with local regulations.

11. Information and Instructions on Disposal and Repairs

Packaging that is no longer required must be disposed of at the local waste disposal facility. The packaging is made of environmentally friendly materials that can be used as secondary raw materials.

The equipment, including accessories and batteries, does not belong in your regular household waste. The EU legislation requires its Member States to collect electrical and electronic equipment and dispose of it separately from other unsorted municipal waste, with the aim of recycling it.

In Germany and many other countries, Sartorius takes care of the return and legally compliant disposal of its electrical and electronic equipment on its own. These products may not be placed with the household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators.

For disposal in Germany and in the other Member States of the European Economic Area (EEA), please contact our service technicians on location or our Service Center in Goettingen, Germany:

Sartorius Stedim Biotech GmbH
Service | FT
August-Spindler-Strasse 11
37079 Goettingen, Germany

Phone +49.551.308.0
Fax +49.551.308.3289

www.sartorius-stedim.com

In countries that are not members of the European Economic Area (EEA) or where no Sartorius affiliates, subsidiaries, dealers or distributors are located, please contact your local authorities or a commercial disposal operator.

Prior to disposal and/or scrapping of the equipment, any batteries should be removed and disposed of in local collection boxes.

Sartorius Stedim Biotech, its affiliates, subsidiaries, dealers and distributors will not take back equipment contaminated with hazardous materials (ABC contamination) – either for repair or disposal. Please refer to the accompanying leaflet | manual or visit our Internet website (www.sartorius-stedim.com) for comprehensive information that includes our service addresses to contact if you plan to send your equipment in for repairs or proper disposal.

12. Appendix

12.1 Detailed Test Description

12.1.1 Diffusion Test (Only valid for the Sartocheck® 4 plus Filter tester)

12. Appendix

When pressure is applied to the upstream side of a wet filter, a diffuse gas stream starts to flow through the filter membrane after a short time.

This gas diffusion can be defined as follows:

$$N = D \cdot L \cdot \frac{\Delta p}{d} \cdot F$$

- N = Molar flow per second
- D = Diffusion coefficient
- L = Solubility coefficient of the test gas in the wetting fluid
- p = Differential pressure available at the membrane
- d = Liquid layer thickness in the membrane
- F = Effective filter area.

The equation above shows that the diffusion also depends on the solubility coefficient "L" and the effective filter area "F," among other factors.

This means that different wetting fluids and test gases will yield different diffusion values on account of the different solubility coefficients, and extremely low diffusion values will be obtained for small filter surfaces (e.g., disc filters).

As a result, these low diffusion values cannot be determined with sufficient accuracy by state-of-the-art measuring technology. This is why the Diffusion Test should not be run on disc filter systems. On the other hand, this test method is ideal for integrity-testing medium- and large-area filter cartridge systems with approximately 1000 cm² and above.

The diffusion of a filter element can be calculated as follows from the pressure drop, test time, upstream system volume and the reference pressure ("p₀" = 1000 mbar), according to the German industrial standard DIN 58356, Part 2:

$$D = \frac{p_1 \cdot V_N}{P_0 \cdot t} \cdot \ln \frac{p_1}{p_1 - \Delta p} \text{ [ml/min]}$$

- D = Diffusion in ml/min or ml/test time
- p₁ = Test pressure in mbar at the beginning of a measurement
- p₂ = Pressure at the end of the measurement in mbar
- Δp = p₁ – p₂ (pressure at the end of a measurement) in mbar
- V_N = Upstream volume of a filter system in ml
- t = Test time in min.
- P₀ = 1000 mbar (or actual measured barometric pressure).

The diffusion calculated is in reference to a constant test pressure and is independent of the system's volume so that a direct correlation with the actual retention capability of a filter is given.

12.1.2 Bubble Point Test (Only valid for the Sartocheck® 4 plus Filter tester)

In wet membranes, the wetting liquid is held inside the porous structure by capillary force. This force increases as the pore size decreases. A specific gas pressure dependent on the pore size of the membrane is necessary to force out the liquid from the pores. This pressure is generally designated as the "Bubble Point," or B.P. for short. It indicates the largest pores of a filter since, liquid is first expelled through them. An overproportional increase in the diffusion rate is characteristic for this process. If the B.P. is below a defined limit, it can no longer be ensured that a filtration process will do the job intended (e.g., remove all bacteria in sterile filtration).

For an ideal capillary system, the Bubble Point can be calculated as follows:

$$\text{B.P.} = \frac{2\sigma \cdot \cos}{r}$$

B.P. = Bubble Point
 σ = Surface tension of the wetting liquid (sigma)
 $\cos$ = Contact angle between the wetting liquid and the filter surface (theta)
 r = Pore radius.

Hence, it follows that the Bubble Point of a filter depends both on the filter material (the contact angle varies from material to material) and, in particular, on the surface tension of the wetting fluid.

Low surface tensions, as are present in organic solvents, detergents and emulsifiers, yield low B.P.s.

If the test pressure required for the diffusion or the Pressure Drop Test is attained, or even if the pressure goes below this level, it is no longer possible to perform the Bubble Point Test under standard conditions. In this case, the test parameters must be changed. This must be kept in mind when the test parameters are entered.

12.1.3 Pressure Drop Test (Only valid for the Sartocheck® 4 plus Filter tester)

The Pressure Drop Test is based on the change in pressure within a filter system connected upstream. This change is caused by the diffusive gas stream flowing through a wet membrane at a given test pressure (see "Diffusion Test"). The pressure drop measured in the process can be directly used as the measure of the filter system's integrity, provided that the upstream volume of the system is constant and known.

At the same diffusion rate, lower pressure drops are yielded for larger upstream volumes, and higher pressure drops for smaller upstream volumes.

Because the result is volume-dependent, this is mainly recommended for leak testing empty housings with a known volume. It can also be derived from the following simplified equation for determining pressure drop based on diffusion, test time, net volume and atmospheric pressure:

$$\frac{D \cdot t \cdot P_0}{V_N} = \Delta p$$

Δp = Pressure drop in mbar within test time "t"

D = Gas diffusion in ml/min

t = Test time in min

P_0 = 1000 mbar or barometric pressure in mbar depending on setting

V_N = Upstream volume of the system in ml.

For this reason, a Pressure Drop Test is recommended mainly to leak-test empty housings and pressure vessels with a known volume.

12.1.4 Plate Holder Test (Only valid for the Sartocheck® 4 plus Bag tester)

The plate holder test and the installed bag test are based on pressure drop measurement. The eventual change in pressure is caused by the diffusive gas stream flowing through an eventual defect in the plastic film of the bag. The pressure drop measures conformity of the bag set-up, provided that the volume of the system is constant and known.

13.1 Network Concept

13. Electronic Data Transfer via Network Connection

The network solution for the Sartochek® 4 plus Filter tester incorporates the TCP-IP and FTP protocol standards, with data being transmitted via the Ethernet standard. FTP uses the passive mode, that exclusively establishes one single TCP connection.

Additional information:

- ▶ Electronic records are transmitted in ASCII (or if selected PDF) format, with an appropriate checksum (CRC).
- ▶ Sartochek® 4 plus is the communication master. All data transfer (import or export) must originate from the Sartochek® 4 plus.
- ▶ Inside the Sartochek® 4 plus, an FTP-client program is running, which actively transmits the electronic records to and from the server PC.
- ▶ Inside the server PC, an FTP-server program is running, which passively receives the electronic records from the Sartochek® 4 plus or provides them to the Sartochek® 4 plus.
- ▶ One server PC (with FTP-server) can serve one or more Sartochek® 4 plus's.
- ▶ One or more local PC's can use the file share.

The following diagram describes the architecture and data flow:


Fig. 13-1 | Diagram showing the architecture and data flow

The following table shows the transmission of data between the Sartochek® 4 plus (SC4 plus) and a PC:

Entity (electronic record)	SC4 plus to PC (export)	PC to SC4 plus (import)
Test Programs	X	X
Test Results	X	X
Log Files	X	
Backup Files	X	X
Configuration Data	X	X

Network Security

The network concept given above is designed for usage in a "closed system"- network on the customer premises. This design concept leads to the following features of the communication setup

- The FTP-data transfer is not encrypted in order to allow the network administrator the customer 's premises to control the network traffic.
- Since the FTP-password and -ID are not encrypted either and could therefore be recorded, these ID-data are system parameters and have no logical link to the user-ID data. Therefore, the user-ID data are still protected.
- To ensure the data integrity of the internal data produced by the Sartochek[®] 4 plus Bag tester, the Sartochek[®] 4 plus Filter tester is the communication master. No remote access is possible on the Sartochek[®] 4 plus system.
- Additionally, a firewall is installed on the Sartochek[®] 4 plus Filter tester to prevent any remote access via network. The firewall allows the Sartochek[®] 4 plus to respond to a network "ping".

13.2 FTP Server Requirements

An FTP server program must be installed on at least one PC to facilitate electronic data transfer. The FTP server program can either be part of your PC's operating system or separately installed software. For the configuration of your FTP server, please refer to the FTP server's manual.

13.2.1 How to Configure Your FTP Server


The following topics should be considered during the configuration of the FTP server:

- Should you create only one or more Sartochek[®] 4 plus Filter tester accounts on your FTP server? This it depends on how you would like to archive your test data. If you would like to save all the data in one folder, you should have only one Sartochek[®] 4 plus account on the FTP server which can be used by all Sartochek[®] 4 plus 's in the network.
- Should it be allowed to overwrite files on the FTP server? In order to avoid the accidental deletion of your data, you should activate the overwrite protection on your FTP server at minimum for the test result folder.
- Different access rights for the different folders? It might be helpful to allow overwriting or deletion for the SC4 plus program folder or SC4 plus config folder on the FTP server if you would like to share this files with several other SC4 plus in the network. Modifications to the files could be carried out on one Sartochek[®] 4 plus and copied afterwards into the FTP folder on the network. Existing older versions will then be overwritten and the newest version will be available on the net work to the other SC4 plus's.

14. PLC Control

14.1 Digital Inputs

There are three digital inputs available (START, PROG and STOP). An impulse at one of these inputs has to be at least 450 msec long for it to be correctly identified.

14.1.1 START Input

An impulse at the START input will start the test program stored in PLC1 or PLC2, depending on the PROG input below (see section 7.13.1 for loading test programs into memory locations PLC1 and PLC2).

If a test unit is currently active, then an impulse at the START input will have no effect.

14.1.2 PROG Input

In order to start PLC2, the PROG input needs to be set to high before the START impulse is initialized. After the test has been completed, the PROG input has to be set to low again.

14.1.3 STOP Input

An impulse at the STOP input has the same effect as pressing the STOP key on the Sartochek® 4 plus Filter tester unit and therefore will stop already started programs on Sartochek® 4 plus Filter tester.

14.2 Digital Output

There are seven digital outlets with the following description and function:

14.2.1 POWER on

This outlet is always high when the PLC-function is activated in the unit.

14.2.2 TEST RUNNING

The outlet is set to high when a START impulse is initiated from the PLC and the additional conditions described above are fulfilled.

14.2.3 ERROR

The ERROR outlet is set to high with every error which can occur during the test run. Additionally, it is high when the temperature control is addressed or the operating pressure in the functional control is too low.

14.2.4 TEST OK

The outlet is high after the test is finished and has passed.

14.2.5 TEST not OK

The outlet is high after the test is finished and has failed.

14.2.6 External Valve 1


The outlet is high if a START impulse is sent by the PLC, and is set back at the end of the measurement phase.


14.2.7 External Valve 2

The outlet is high if a START impulse is sent by the PLC, and is set back at the end of the measurement phase.


All outlets are chargeable with 0.3 A. All outlets except the "POWER ON" display are reset at high with a STOP impulse. After this STOP impulse, the unit is again ready for the next test.

14.3 Pin Layout



PLC Connection

Simple Relay Connection

12 pin Connector, male
Solder side

14.5 Outlets

Name	Function	Outlet Status						
		D	E	F	G	H	L	M
Outlet D	Unit is ready to use (stand by)	1	0	0	0	0	0	0
Outlet E	Test running	1	1	0	0	0	1 ¹	1 ²
Outlet F	Error display	1	1	1	0	0	0	0
Outlet G	Test passed (OK)	1	1	0	1	0	0	0
Outlet H	Test failed (not OK)	1	1	0	0	1	0	0
Outlet L	Ext. valve 1	see note 1						
Outlet M	Ext. valve 2	see note 2						

Note 1: Outlet L not activated during a test until the end of stabilization time

Note 2: Outlet M not activated during a test until the end of measuring time

Plug Layout

Plug	Function
Pin A	START
Pin B	STOP
Pin K	PROG
Pin C	GND
Pin D	POWER ON
Pin E	TEST RUNNING
Pin F	ERROR
Pin G	TEST OK
Pin H	TEST NOT OK
Pin J	+24 V DC
Pin L	EXT. VALVE 1
Pin M	EXT. VALVE 2

15. Index

Ambient conditions and conditions at the place of installation	9, 16	Hardware configuration	53f, 57
Appendix	7, 31, 35ff, 81	Hazard(s)	10, 14
Auto Log Off	60f	Housing	6, 11, 17, 21, 23f, 26, 29, 31, 33f, 36, 38, 51, 52, 70, 76ff
Backup all test programs on to SD card (F3)	42	If the delivery has other defects	16
Bubble point test	30f, 35f, 58, 59, 76, 82	If the shipment is not complete	15
Calibration	24, 29, 53, 62, 63f, 65, 71f	Incoming inspection	15
Calibration adjustment (only with multi test mode)	62f, 65	Insert paper roll and ribbon cartridge	21
Changing the fine-wire fuses	68	Installation	5, 7, 13ff, 39
Cleaning	12f, 14, 24, 50f, 68, 75	Installation instructions	7, 16
Cleaning (F3)	24, 50	Installation of replacement or wear parts	14
Cleaning drying	24, 60	Installed safety equipment	11
Complaints	15	Interface port for data output	8
Connection data	7, 28	Interim storage	16
Connections	5, 10, 17f, 20, 28, 32, 50ff	Internal pressure sensor	18f, 24, 52f, 62f
Connections to the internal pressure sensor	18	Left and right panels	17f
Connections to the external pressure sensor	20	Line disconnecter	11
Contents	2ff	Load test program from SD card (F2)	42
Copyright	2, 94	Load test program from internal memory (F1)	42
Diffusion test	30ff, 36ff, 76f, 81, 83	Loss of compressed air	29
Diffusion test and bubble	30	Maintenance	7, 11ff, 68, 69
Dimensions and weight	7	Maintenance servicing	68
Disposal	12, 80	Manage users	66
Drying (F4)	52	Manufacturer's ID label	10
Emergencies	80	Measuring accuracy	8
Equipment supplied	15f	Measuring ranges	8
External pressure sensor	18, 20f, 32, 39, 52f, 62, 71	More detailed test description	81
Front panel and back panel	17	Multipoint diffusion test	30, 38, 39, 77
FTP Server	43, 55f, 56, 85	Network connection	8, 17, 55f, 74, 84
Function	11, 19, 23ff, 43, 45, 47, 48ff, 57, 61f, 67, 74ff	Notes	10, 13, 16
Function test	19, 24, 49, 50	Notes explanations	10
Function test MU (F1)	49	Operating conditions	8, 13
Function test TU (F2), selection test unit	50	Operating elements	26
General data	9	Operating personnel and service and Maintenance personnel	11, 13
General functional description	23	Operation	7, 9, 13, 18, 26ff
General warnings	13, 14	Operator's obligations	12
Grounding conductor connection	10	Overview	6ff, 23, 26, 58
		Overview and proper use	6ff

Packaging for return shipment	16	Technical specifications	7, 16, 28
Panels	17f	Test media	7
Plug connections	10	Test methods	9, 11, 30
Pneumatics	7	Test parameters	29ff, 37ff, 70, 82
Power failure	24, 29f	Test programs	23, 42, 43, 45f, 58, 67, 74
Pressure drop test	30, 38, 76f	Test sequence	34, 37, 40f, 68
Pressure gauge function	24, 49, 52	Tests	6, 12, 23f, 27ff, 39, 49, 60
Product specific data	8	Touch screen	6, 8, 24, 26ff, 30f, 43, 48, 49, 71
Program database (F4 – Main Menu)	42	Trained workers	11, 69
Programming tests (F2 – Main Menu)	30, 31	Translation	2
Proper use	6ff	Transporting and unpacking the equipment	16
Protective equipment	11	Troubleshooting	13, 69
Receipt of delivery	15	User interrupt	24
Results database (F5 – Main Menu)	44	Visualize (F3 – Main Menu)	47
Safety	10ff, 16, 24, 69f	Water flow test	37, 72
Safety functions	24	Water intrusion test	29f, 37, 70, 72, 77
Safety precautions	11ff, 69		
Safety tests and inspections	12		
Selectable languages	8		
Self test for broken conductors and short-circuiting	24		
Serial number	11, 25, 48, 57, 63, 69		
Service (F7 – Main Menu)	53		
Service Menu	27, 32, 51f, 58ff, 67f, 70ff		
Shutdown procedure	14, 68f		
Software configuration	51, 53, 58ff		
Software update	25		
Special functions (F6 – Main Menu)	46		
Switching off the unit	11, 29, 30		
Switching off the unit in the event of a power failure	29		
Switching the unit ON	28		
Switching the unit ON and OFF	28		

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Phone: +49.551.308.0
Fax: +49.551.308.3289
www.sartorius.com

The information and figures contained in these instructions correspond to the version date specified below.

Sartorius reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

Copyright notice:

This instruction manual, including all of its components, is protected by copyright. Any use beyond the limits of the copyright law is not permitted without our approval. This applies in particular to reprinting, translation and editing irrespective of the type of media used.

© Sartorius Germany

Last updated:
09 | 2015